



# 371, 311 & HB80 Medisafe QCCP Vaccine Fridges

## Description

Fridges specifically designed for vaccine storage. Temperature control system pre-configured to hold temperature from 2°C to 8°C (approved QCCP specification). Easily download temperature history to any computer via the data logger. Available in three sizes.

## Features



- Electronic temperature control and monitoring system
- Audible auto resetting alarms
- External Hi Lo Data Logger
- Computer Data Logger
- Auto dialer (option)



COOL SOLUTIONS

## Dimensions

Model	H (mm)	W (mm)	D (mm)	Litres
371	1930	595	595	381
311	1630	595	595	306
HB80	855	500	560	80



## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Q 4077 PO Box 965 Sumner Park BC Q 4074  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au www.medsurg.com.au Fax 07 3859 2995

