

All offers expire
31/08/06

MEDI+Saver

All prices plus GST unless
otherwise indicated

July/August SPECIALS

NEW PRODUCT

Introducing the new CARDIOLINE® range of ECG's.

Features of the new CARDIOLINE® range include

- Automatic, manual and pre-programmed recording mode - LCD display - Alphanumeric keyboard - Dual power supply: mains and rechargeable internal batteries - Memory of the last ECG recording and printing of additional copies - Safe and easy wireless connection, thanks to the IR (Infrared) digital interface - Time and date indication

The CARDIOLINE® ar600adv

A portable electrocardiograph and much more
The CARDIOLINE® ar1200view and ar2100view
Easy view... clear recordings



The CARDIOLINE® ar600adv represents the evolution of a portable instrument able to meet the growing demands placed on ECG diagnostic systems.

(Multi-channel print format on 60mm z-fold paper: 1,2, 2+R, 3 channels)

EZ-AR600ADV *From \$3.20 per day (2 years)

\$2,091.00 *GST FREE



The CARDIOLINE® ar1200view is a multi-channel electrocardiograph, providing solutions that simplify the activity performed by specialists and general practitioners alike.

(Multi-channel print format on 120mm z-fold paper: 3, 3+R, 6 channels)

EZ-AR1200VIEW *From \$3.95 per day (2 years)

\$2,579.00 *GST FREE



The CARDIOLINE® ar2100view, in short, is an A4 format electrocardiograph providing solutions that will help improve the quality of diagnosis performed by specialists and general practitioners alike.

(Multi-channel print format on A4 paper: 3, 6x1, 6x2, "Full Page" (3x4+R) and (3x4+3R) 12 channels)

EZ-AR2100VIEW *From \$6.07 per day (2 years)

\$3,921.00 *GST FREE



Sunnex® Surgical and Examination Lighting

The SoloSpot is ideal for applications where a focused, spot light is required. The specially designed lamp head offers cool light and a focused, intense spot. Ideal for OB-GYN, dental and surgical applications.

NEW TO MEDICAL & SURGICAL



Features include

- Unique double filtered lens and colour correction filter with lampshade extension, helps keep light cool on skin and illuminate with natural skin tones to improve tissue clarity
- Excellent manoeuvrability allows for easy adjustment of light
- Intense 20W Halogen light source (average 2,000 hours)
- 3.5 inch diameter spot pattern @ 20 inches

Sunnex SoloSpot™ Halogen Light - Wall mounted (White)

PL-HM2005-37EXT

\$1,010.00

Sunnex SoloSpot™ Halogen Light - on Mobile Base (White)

PL-HM2006-38EXT

\$1,045.00



The Sunnex® HM Series Multipurpose Examination Light

offers intense white light in a compact, flexible design. Features include

- Compact, vented shade minimises heat output
- Stay cool finger guard for easy adjustment of light
- Swivel joint lamp head allows 359° axial rotation and 90° angular adjustment

Sunnex HM Series Halogen Light - Wall mounted (White)

PL-HM1000-37W

\$308.75

Sunnex HM Series Halogen Light - on Mobile Base (White)

PL-HM1006-37W

\$395.00

Medilite Wall Mounted X-Ray Viewers



Manufactured from high quality aluminium with a modern and attractive rounded styling, Medilite x-ray viewers are durable and lightweight and are designed to meet the high standards of the medical profession.

- Unique top film-grip holder design ensures a light yet positive grip on films and is distortion free
- Durable white anodised aluminium
- Internal baffle system for even light distribution

Single bay wall mounted x-ray viewer, 370mmW x 500mmH x 110mmD

MT-001W

\$129.15

Double bay wall mounted x-ray viewer, 720mmW x 500mmH x 110mmD

MT-002W

\$258.30

NEW DRESSINGS TROLLEYS

Exclusive to Medical and Surgical



Stainless Steel Dressing Trolley 500mm x 500mm x 900mm. With top and bottom rail and two lockable castors

FANTASTIC BUY!

J1-1-001

\$365.00

Check out our range of vaccine refrigerators

* Please note that freight charges will apply

The G135L Vaccine Fridge



Thermoline Scientific range of vaccine fridges



NEW TO MEDICAL & SURGICAL

Thermoline vaccine refrigerators will suit situations where pharmacy/vaccine products are required to be stored safely. Temperature is precisely controlled electronically at 5°C in ambients up to 32°C while an independent digital indicator of temperature contains the Hi (+8°C) / Lo (+2°C) temperature alarms emitting visual, audible and remote signals. All Thermoline Economy Pharmacy/Vaccine Refrigerators incorporate the following features:

- 595mm W x 600mm D x 830mm H
- 135 Litre capacity
- Temperature range of +2°C to +8°C @ 30°C ambient temp + 60% humidity
- Adjustable shelves
- Audible alarm system with additional external alarm connection
- Ventilated cooling system
- Electronic temperature control + display
- Door lock
- Automatic defrost
- Maintenance free condenser
- * Fully compliant with Quality Care Pharmacy Program (QCPP) Cold Chain Standards
- * 3 year warranty

- Internal forced air flow
- Separate digital temperature display with alarms and temperature controller
- Max/min temperature logging
- Door keylock
- Adjustable/removable shelves
- Remote monitoring outlet
- Glass doors

140 Litre, 890mm x 600mm x 610mm, 3 shelves

VZ-TEPR140-1-GD

\$2,925.00

250 Litre, 1690mm x 590mm x 620mm, 4 shelves

VZ-TEPR250-1-GD

\$3,362.00

400 Litre, 2030mm x 700mm x 610mm, 5 shelves

VZ-TEPR400-1-GD

\$3,550.00

\$1,325.00



MEDICAL AND SURGICAL

Requisites Pty Ltd

Established in 1922



Check us out Online!
Visit our website to view our online catalogue. www.medsurg.com.au

ABN 70 009 659 394

33 Fulcrum Street, Richlands Q 4077

Phone 07 3217 0966 Fax 07 3375 1744 Toll Free 1800 801 226

Email: medical@medsurg.com.au

www.medsurg.com.au

Rep Talk



Medfin Australia

Medfin are specialists in providing tax effective, low cost finance structures to the medical, dental and allied health professionals. Medfin's range of financial products are specifically designed to simplify the process and take the hassle out of arranging your finance so you can spend more time in your practice and doing the things you enjoy.

Medfin has established itself as Australia's leading financier to the medical, dental and allied health professions. Medfin's reputation relies on providing exceptional service and tailored products, efficiently and professionally, as a viable alternative to traditional banking.

Medfin will fund the assets for your practice including:

- * Equipment - medical/dental
- * Practice purchases
- * Fit out
- * Lines of credit
- * Computer equipment and software

Medfin have a flexible approach which allows you to choose the optimal structure to suit your needs. Medfin Finance: Fast approvals, no forms to complete - It couldn't be easier!

Kimberley Thiele
Territory Manager



The Pony FX Desktop Spirometer

The Pony FX designed by Cosmed for lung function screening, includes all of the experience, latest requirements and customer wishes collected in more than 20 years. * Medical Director Compatible



Key Features include

- Screening spirometry (FVC - SVC - MVV)
- Large graphic colour display
- High speed thermal printer (ATS compatible)
- Battery and mains operating

Pony Fx Spirometer

SZ-42FX *From \$4.05 per day (2 years)

\$2,634.50



Accessories for the Pony FX

Paper for the Pony FX - 112mm roll

SZ-43

\$12.94

Bacterial Filter for Pony Spirometers (Disposable)

Packets of 50

SZ-PONYFILTER

\$126.75

Disposable mouthpieces for Pony FX Spirometer

Adult - Box of 200

SZ-42D

\$42.00

Disposable mouthpieces for Pony FX Spirometer

Child - Box of 200

SZ-42DD

\$38.40

Printer ribbon for Pony FX

SZ-42E

\$35.40

Soniclean

A new generation of ultrasonic cleaners

There's no faster, safer, quieter or easier way to thoroughly clean items than by using a Soniclean ultrasonic cleaner. Soniclean provides environmentally safe, highly efficient, splash-free cleaning power with an absolute minimum of energy consumption.

- Completely sealed units - meets IP 66 and AS2773
 - Keypad touch control
 - Over temperature safety feature
 - Corrosion resistant powder-coated metal case
- * Please note that "T" indicates Timer and "HT" indicates Heater and Timer



2 Litre Soniclean - overall size 183mm x 145mm x 172mm high

SZ-120T

\$736.00

Perforated S/S Tray/Basket for 2L

SZ-SA120

\$78.00

S/S Lid/Cover for 2L

SZ-SC120

\$16.90

3 Litre Soniclean - overall size 273mm x 145mm x 172mm high

SZ-160T

\$854.00

Perforated S/S Tray/Basket for 3L

SZ-SA160

\$90.00

S/S Lid/Cover for 3L

SZ-SC160

\$23.00

6 Litre Soniclean - overall size 333mm x 200mm x 188mm high

SZ-250T

\$1,244.00

Perforated S/S Tray/Basket for 6L

SZ-SA250

\$108.00

S/S Lid/Cover for 6L

SZ-SC250

\$30.00

10 Litre Soniclean - overall size 333mm x 200mm x 275mm high

SZ-500HT

\$2,269.00

Perforated S/S Tray/Basket for 10L

SZ-SA500

\$126.00

S/S Lid/Cover for 10L

SZ-SC500

\$36.00

After some scales?

Proport® will help you weigh up your options!

Electronic Baby Scale



- 20kg or 44lb capacity - Large power saving LCD readout - Lock in weight feature
- a stable weight can be obtained even if the baby is moving - Removable tray for easy storage or weighing of a small child
- Requires only one 9 volt alkaline battery (not included) - 12 month warranty

PRO-1600

\$162.00

Precision Bathroom Scale



- 136kg or 21st 6lb capacity
- Displays weight in both kilograms or stones/pounds - Large easy-to-read 25mm LCD - Requires only one 9 volt alkaline battery (not included) - 12 month warranty

PRO-1650

\$44.75

Vogue Bathroom Scale



- 150kg or 23st capacity
- Displays weight in both kilograms or stones/pounds - Large easy-to-read 32mm LCD with backlight (easily visible under bright or dim lighting) - Stainless steel platform - Requires only one 9 volt alkaline battery (not included)
- 12 month warranty

PRO-1750

\$59.75

Big Foot Bathroom Scale



- 150kg or 23st capacity - Displays weight in both kilograms or stones/pounds
- Very large foot platform with easy-to-read speedometer type dial
- 5 year warranty

PRO-2800

\$37.50

VALUE!

An offer that won't leave you breathless!

The Bestneb and Liberty® nebuliser pumps

The Bestneb Nebuliser Pump has been designed with children in mind to make respiratory therapy a less frightening experience. Along with a super quiet operation, the Bestneb Nebuliser Pump should be your first choice when considering a nebulisation unit for the home or surgery!



Bestneb Nebuliser Pump

LR-PN

\$119.90 *GST FREE

The Liberty® Nebuliser Pump is the result of 20 years of extensive research and development to provide the patient with the most cost effective, versatile and reliable pump on the market.



Liberty® Nebuliser Pump

LR-PNL

\$78.00 *GST FREE

FASTER, SAFER, QUICKER!

Dalcross IV Stands with mobile base

- Stable 5-star solid steel base, finished in hammertone gunmetal grey powder coat
- Chrome plated stem with keywheel for height and tension adjustment
- Black 50mm twin wheel castors



Mobile IV Stand with stainless steel pole

DL-D19DS/S
\$206.50

Mobile IV Stand with chrome pole

DL-D19DS
\$187.50

Medsure Latex Non-sterile Examination Gloves

Medsure Powdered Gloves are manufactured with a beaded cuff and contain a thin coating of cornstarch powder for easy handling and to prevent sticking. Medsure Powdered Gloves have a low powder and protein content for minimalised sensitisation, and are of a high strength and elasticity, to resist tearing.



Medsure Powder Free Gloves undergo extensive washing and treating by a chlorination process. This results in the elimination of most latex protein.

Available in Boxes of 100

MS-1-01 (Small)
MS-1-02 (Medium)
MS-1-03 (Large)

\$4.85 1 - 9 Boxes
\$4.35 10+ Boxes

Available in Boxes of 100

MS-1-04 (Small)
MS-1-05 (Medium)
MS-1-06 (Large)

\$6.00 1 - 9 Boxes
\$5.80 10+ Boxes

Also available - 1 and 3 tiered glove dispensers

MD-GD100 (1 tiered)

\$5.55

MD-GD101 (3 tiered)

\$21.15

IV Solutions



Sodium Chloride 0.9% for I.V. - 2x100ml bag
TV-AHB1364

\$8.10

Sodium Chloride 0.9% for I.V. - 500ml bag
TV-AHB1323

\$2.65

Sodium Chloride 0.9% for I.V. - 1000ml bag
TV-AHB1324

\$3.85

0.18% Sodium Chloride and 4% Glucose for I.V. - 500ml bag
TV-AHB1253

\$2.68

0.18% Sodium Chloride and 4% Glucose for I.V. - 1000ml bag
TV-AHB1254

\$3.85

5% Glucose for I.V.- 500ml bag
TV-AHB0063

\$2.65

5% Glucose for I.V.- 1000ml bag
TV-AHB0064

\$3.85

Ringers solution - 1000ml I.V. bag
TV-AHB2304

\$5.85

Compound Sodium Lactate (Hartmann's solution) for I.V. - 500ml bag
TV-AHB2323

\$2.65

Compound Sodium Lactate (Hartmann's solution) for I.V. - 1000ml bag
TV-AHB2324

\$3.85

FOLIODRAPE® A cost effective drape sheet from Paul Hartmann

45cm x 75cm, available in boxes of 65



Foliodrape® plain drape

HA-250490/1

\$47.60

Foliodrape® fenestrated drape (non-adhesive) two layer
HA-250493/1

\$58.72

Foliodrape® fenestrated drape (adhesive) four layer
HA-250497/1

\$122.40

Introducing the Breath-A-Tech spacing device and child mask

- Hospital grade quality (will withstand autoclaving)
- Can be used with standard puffers in the treatment of asthma
- The mask attaches easily to the spacer
- Makes the administration of a child's medication easier and more effective

* Please note that the mask and spacer are sold separately

Breath-A-Tech hospital grade spacer

BZ-BAT06

\$20.15 * GST FREE

Silicone child mask for use with Breath-A-Tech spacer

BZ-BAT07

\$11.70 * GST FREE

NEW TO
MEDICAL & SURGICAL



THE BIO-CAN® RANGE OF SHARPS CONTAINERS

For the cost effective, safe containment of sharps, biohazardous and clinical waste

* Manufactured to conform to AS-4031/92



125 ml Yellow round BIO-CAN® with screw cap
PI-742007

\$0.68

250 ml Yellow round BIO-CAN® with screw cap
PI-742008

\$0.93

600 ml Yellow round BIO-CAN® with screw cap
PI-742001

\$2.10

1.4 Litre Yellow square BIO-CAN® with funnel top and press down lid
PI-1-01

\$3.69

3.2 Litre Yellow oval BIO-CAN® with push cap
PI-7420022

\$7.39

5 Litre Yellow oval BIO-CAN® with hand guard and funnel top
PI-7420023

\$7.90

8 Litre Yellow oval BIO-CAN® with hand guard and funnel top
PI-7420024

\$9.30



Surgical Paper Tape - Hypo-allergenic, white
12.5mm x 9m - Box of 24
EG-3005109005

25mm x 9m - Box of 12
EG-3005109010

\$7.95



Underpads - 5 ply, 68.5cm x 42.5cm, Box of 250 (5 packets of 50)
EG-4818402030

\$27.50

1-2 Boxes * GST FREE

\$26.50 3+ Boxes



Wooden Tongue Depressors, 15cm - Box of 100
EG-4421909090

\$1.00

New style ECG Universal multi clip adaptors

For tab and snap electrodes, these new style clips will fit 3mm to 4mm pin leads

Packets of 10
3M-UNI

\$21.07



Dry Comfort Plus Bed Pads
860mm x 910mm with wings and without waterproof backing (White)
MB-820565W

\$27.00 * GST FREE

860mm x 910mm with wings and with waterproof backing (Green)
MB-820565WB

\$30.10 * GST FREE



Plastic Disposable Aprons
71 x 117cm, packs of 100
MPM-AP001/100

\$10.15

Plastic Disposable Aprons
71 x 117cm, 10 packs of 100
MPM-AP001

\$85.40

Plastic Disposable Aprons
70 x 115cm, Tear off, 2 packs of 50's to hang on wall
MPM-CA-APRON-1

\$9.15

Blue Plastic Overshoes
packs of 100
AC-300050/100

\$5.50

Blue Plastic Overshoes
10 packs of 100
AC-300050

\$36.00



"Medical and Surgical Requisites Quality products at a quality price!"

At MSR we understand market trends. In response to this we source products from many different countries, ensuring that we have the best selection to provide you with the most suitable and cost effective products. This means that we have become one of Australia's largest importers of exclusive medical and nursing products.



Medical and Surgical Requisites Pty Ltd
"We share the responsibility for your patient care"

Gary Monaghan
Territory Manager

Save time and money with MEDI+Order®

MEDI+Order is an interactive ordering service designed to eliminate ordering mistakes, while saving you time and money.

We simply send you a list of the products your practice or nursing home has ordered in the past six months. You mark the products you require along with the quantities you need and send it back. Your MEDI+Order will be emailed or faxed to you on whatever date you choose. You simply look down the list, mark the products you want with the quantity and email, fax or post it back to us.

CALL 1800 801226

No Account..... No Problem!

Introducing...



As Managing Director, I am proud of Medical and Surgical Requisite's established reputation since our commencement in 1922.

Geoffrey Bailey
Managing Director



The MSR / Paul Hartmann Continence Care Programme is the ultimate solution, specifically designed around the needs of your organisation.

Maureen Cadwallander
Manager, MSR Continence Care Division



Our dedicated sales team is focused on new products and procedures to meet market demands.

Kathy Cubis
Sales Manager



MEDI+Plus is focused on providing your surgery with unique support and the widest possible range of products. We regularly inform you of new services, products and cost saving benefits.

Lorraine Brown
MEDI+Plus Consultant



The MEDI+Gate stock control system has been designed to simplify the ordering of continence aids for your nursing facility.

Courtney Barrett
Customer Support and Systems Coordinator

**We will have your new account
open on the same day. We also
accept all major credit cards.
Simply place your order:**



By Email:

For simple and easy processing of your order at the click of a button, email your orders to:
medical@medsurg.com.au



By Fax:

Send your order by fax.

Fax No: (07) 3375 1744



By Phone:

Just phone and place your order or discuss your needs. We will assist you to directly place your order over the phone.

Free Call: 1800 801 226



In Person:

Discuss your needs directly with your territory manager, either by phone or in person. Our representative will be pleased to call as required.

Telephone: 07 3217 0966



By Mail:

We still receive many orders by post. This can be very convenient especially if you need to send an order on your pre-printed form.

Medical and Surgical Requisites Pty Ltd
PO Box 166
Richlands Q 4077 Australia

**Credit cards do not apply to account
customers. Cash sales only.**



**MEDICAL
AND
SURGICAL**

Requisites Pty Ltd

Established in 1922



Check us out Online!
Visit our website to view our online
catalogue. www.medsurg.com.au

ABN 70 009 659 394

33 Fulcrum Street, Richlands Q 4077

Phone 07 3217 0966 Fax 07 3375 1744 Toll Free 1800 801 226

Email: medical@medsurg.com.au www.medsurg.com.au