

Introducing the next generation of a clinical icon.



### Open bell versatility.

The paediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim.

## Introducing the 3M™ Littmann® Classic III™ Stethoscope

### Right out of the box, it lives up to its name.

The Littmann Classic III stethoscope is the latest version of the stethoscope that helps millions of medical professionals achieve their best. It brings new design, materials, and technology to the Littmann stethoscope series trusted by clinicians around the world. The new Classic III stethoscope offers you:



**High acoustic sensitivity** for performing general physical assessments.



**Tunable diaphragms** on both the adult and paediatric sides of the chestpiece.



**Single-piece diaphragm** that's easy to attach and easy to clean because the surface is smooth without crevices.



**Next-generation tubing** with improved resistance to skin oils and alcohol for longer tubing life. Less likely to pick up stains.

Using a Classic III stethoscope demonstrates what's important to you as a medical professional: quality, achievement, and providing the best possible care. Every patient. Every day.



Established 1922

### Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
Phone 1300medsurg (1300 633 787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995



Here's what makes the

# 3M™ Littmann® Classic III™ Stethoscope stand out.

## Features

**Tunable technology** is a 3M invention that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece.

**Soft-sealing eartips** give you a comfortable fit and an excellent acoustic seal; they snap tightly onto the eartubes for safety.

**Open bell stays clear** of dirt and debris by keeping it covered with the small diaphragm.

**Adjustable headset** lets you easily set the tension for head size and comfort just by squeezing or pulling apart the eartubes.

**No natural rubber latex or phthalate plasticisers** are used in the tubing or any other component—to help protect human health and the environment.

### Classic III stethoscope specifications

**Length** 69cm / 27.2 inches

**Weight** 150g

**Acoustic performance** 7

#### CHESTPIECE

**Construction** Two-sided, adult and paediatric

**Material** Stainless steel

**Diaphragm type** Tunable, single-piece

**Diaphragm diameter** Adult: 4.3cm / 1.7 in  
Paediatric: 3.3cm / 1.3 in

#### TUBING

**Construction** Single lumen (single acoustic channel within tube)

**Material** Not made with natural rubber latex or phthalate plasticisers

#### HEADSET

**Material** Wide diameter aerospace alloy / anodized aluminium

### Available tubing colours and finish combinations

TUBING COLOUR	STANDARD FINISH	SPECIAL FINISH	TUBING COLOUR	STANDARD FINISH	SPECIAL FINISH
Black	5620	Black 5803, Smoke 5811	Grey	5621	
Navy Blue	5622		Caribbean Blue	5623	Rainbow 5807
Raspberry	5626	Rainbow 5806	Hunter Green	5624	
Burgundy	5627		Orange	5629	
Chocolate		Copper 5809	Ceal Blue	5630	
Lavender	5832		Peal Pink	5633	
Turquoise	5835		Lime Green	5829	
Lemon-Lime	5839		Plum	5831	
Emerald	5840		Dark Olive Green		Smoke 5812
Rose Pink	5639				

## Colours and Finishes

We periodically introduce new tubing colours and chestpiece finish combinations for the Classic III stethoscope. Visit either [www.littmann.com.au](http://www.littmann.com.au) or [www.littmann.co.nz](http://www.littmann.co.nz) for the latest product information.

*Due to limitations of the printing process, the colors shown are approximations of the actual tubing colors.*



### Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
Phone 1300medsurg (1300 633 787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995

