

A100D Portable Fetal Doppler



CE 0123

*Color OLED Fetal Doppler!
Algorithm inside, More accurate!*



A100D Portable Fetal Doppler is a high performance fetal heart beat detector with color OLED displaying FHR value, which can meet routine examination requirements of obstetricians. It is applicable for the private obstetrician's office, antepartum clinic, obstetrical department of hospital and remote monitoring situations, or using at home.

- Lightweight and Easy-To-Use.
- Color OLED display FHR value, FHR intensity, speaker volume, battery capacity indicator, working mode, frequency of the probe and FHR synchronization symbol.
- Accurate FHR detection. It has FHR mode, average FHR mode, manual count FHR mode and count-down FHR mode.
- Detection and Indication of probe off.
- Improper FHR alarm value adjustable. The blinking FHR digit indicates improper FHR value.
- Turn to power saving mode automatically without signal within 30 seconds.
- Auto shut off without signal within 2 minutes.
- Earphone port supports FHR sound output to earphone or recorder.
- Standard 9V alkaline battery is available for power supply, which can support 300 minutes continuous work.

Features

Specifications

FHR

Ultrasound Frequency: 2.0MHz
Intensity : 5mW/cm² or less
FHR Range: 60-210bpm
FHR Range for manual count mode: 30-250bpm
FHR Accuracy: ± 1 bpm

Display

Color OLED display FHR value, FHR intensity, speaker volume, battery capacity Indicator, working mode, frequency of the probe and FHR synchronization symbol.

Standard Accessory

2.0MHz Doppler Probe 1 pc
Operation Manual 1 pc
Charger 1 pc
Rechargeable 300mAh Ni-MH battery (9V) 1 pc

Optional Accessory

Ultrasound Gel 1 pc
9V Alkaline Battery 1 pc

Dimension: 142mm(L) x 112mm(W) x 42mm(H)

Outer Box Dimension: 210mm (L) x 180mm (W) x 60mm (H)

Net Weight: 0.3 kg

Gross Weight: 0.47 kg

Outer Box Dimension (20 boxes):

630mm(L) x 220mm(W) x 370mm(H)



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

