

Bovie®

SPECIALIST|PRO

Electrosurgical Generator



The Bovie® Specialist Pro is the electrosurgery solution for in demand physician offices and specialty facilities. Designed to operate like an operating room generator in a compact lower watt package. The multipurpose unit is ideal for monopolar and bipolar office surgical procedures.

The Specialist|PRO offers 5 common energy outputs controlled by the new sleek, push-button design, with up to 120 watts of power. The modalities include Cut, Blend, Coagulation, Fulguration and Standard Bipolar.

Industry-standard connections allow for use of a variety of monopolar and bipolar accessories providing energy versatility across a variety of specialties such as OBGYN, Urology, Plastic Surgery, Dermatology, ER and more.

Bovie NEM™ neutral electrode monitoring increases patient safety by ensuring proper pad connectivity during procedures.

The Bovie Specialist|PRO is backed by the unmatched **4-year manufacturer warranty**, offered only by Bovie.

Efficiency, Versatility, and Safety ...
Designed benefits for you and your patients

DALCROSS®
Medical Equipment



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



BOVIE® SPECIALIST|PRO ELECTROSURGICAL GENERATOR



STANDARD EQUIPMENT

- A1250S electrocautery generator,
 - ESP1 sterile disposable pencil
 - ES01 blade, ES02 needle, ES20 ball (1 each) sterile disposable electrodes
 - ESRE (5) split disposable return electrode
 - A1252C (2) reusable grounding cords
 - 110 VAC Hospital-grade power cord*
 - 220 VAC Hospital-grade power cord*
 - User's guide (MC-55-237-001)*
- * Not Pictured



ESP1



ES01, ES02, ES20



ESRE



A1252C

COMPLEMENTARY PRODUCTS †



SE02
Smoke Evacuator †



BV-1253B
Monopolar Footswitch †



A1255A
Footswitch Adapter †




Disposable LLETZ/LEEP
Loop and Ball Electrodes †



Split Disposable Return
Electrode with cord †



Disposable SuperCut™
Needles Electrodes †



Disposable Resistick II™
Electrodes †



Disposable Dermal Loops †

† Sold Separately

PRODUCT SPECIFICATIONS

DIMENSIONS Height:
15.2cm
Width: 26cm
Depth: 30.5cm
Weight: < 4Kg

OUTPUT POWER
Cut: 120 watts
Blend: 90 watts (60% Cut, 40% Coag)
Coagulation: 80 watts
Fulguration: 40 watts
Bipolar: 30 watts

OUTPUT FREQUENCY
350 – 600 kHz
LINE VOLTAGE
A1250S: 100-240 VAC

LINE FREQUENCY
50 – 60 Hz
WARRANTY
Four Years



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsur (1300 633 787) Email medical@medsur.com.au Fax 07 3859 2995

