



The perfect scale for daily use at Hospitals, Clinics and GPs offices. Strong and sturdy construction with a total capacity of 250 kg and a fine graduation of 100 g. The turning display allows the reading from both sides of the display. Mains and battery operated for maximum flexibility.

Technical specification:

Capacity:	250 kg	Function:	Mother and child weighing, turning display,
Graduation:	100 g		auto switch off
Dimensions:	400 x 350 x 44 mm	Power supply:	4 x AA batteries life time approx. 4000 cycles - mains optionally.
Weight:	3,8 kg	Levelling feet	
Main all off switch		Tap on start function	

ADE (GmbH & Co) - Hammer Steindamm 27-29 - 22089 Hamburg - Germany - www.ade-germany.com



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
 Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

