



The perfect scale for daily use at Hospitals, Clinics and GPs offices. Strong and sturdy construction with a total capacity of 250 kg and a fine graduation of 50 g. The turning display allows the reading from both sides of the display. Mains and battery operated for maximum flexibility. A BMI function is integrated as standard function.

**Technical specification:**

Capacity:	250 kg	Function:	BMI calculation, mother and child weighing, turning display, hold & auto-hold, auto-switch off
Graduation:	50g<50kg>100 g	Power supply:	4 x AA batteries life time approx. 4000 cycles and mains
Dimensions:	400 x 350 x 44 mm	Levelling feet	
Weight:	3,8 kg	Tap on start function	
Main all off switch			



Established 1922

**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

