



Quick to use portable height measuring device. The inbuilt spirit level will help to place the height measure correctly on the head of the patient. The measuring is done in seconds via an ultrasound signal and will be shown on the integrated LCD display. Foldable for easy transport.

Technical specification:

Measuring range	100-1000 mm
Graduation	1 mm
Dimensions	1260 x 280 x 90 mm
Weight approx.	585 g



ADE (GmbH & Co) - Hammer Steindamm 27-29 - 22089 Hamburg - Germany - www.ade-germany.com



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
 Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

