

# Heli-Checks

**A hollow load challenge test for pre-vacuum steam sterilizers.**

## Hollow Load Testing:

- Meets EN 867-5 requirements
- A good vacuum for steam penetration at proper time and temperature will change indicator colour from purple to green
- Reusable device for 250 uses
- Simple to load and remove indicator
- Low cost per test
- Package contains one device with resealable



bag of 250 high-precision sterilant specific indicators

- 3-year shelf life
- Helix devices are not useful for testing of porous loads and therefore should not be used as a Bowie-Dick test



**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900  
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077  
Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au)

[www.medsurg.com.au](http://www.medsurg.com.au)  
Fax 07 3859 2995

