

AQ Medical Vaccine Fridges
are a **robust, cost effective**
solution for all of your **vaccine**
refrigeration needs.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



AQ Medical Vaccine Fridges are a **robust, cost effective** solution for all of your **vaccine refrigeration** needs.

Features

- Self closing door to ensure temp never breaches vaccine temps
- Temperature keeps from 2°C - 8°C accurately and evenly
- High-low temp audible and visual alarm
- Suitable for small vaccine bottles
- Factory installed lock ensures security
- Electronic temperature controller and display
- 24 months warranty with national service support



AQM 110



Dimensions (HxWxD)	820 x 595 x 595 mm
Capacity	145L
Temp Range	2°C - 8°C
No. of Shelves	4
Warranty (parts and labour)	24 Months
Power Supply	10amp
Feet Included	Feet at front, rollers at back



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995





AQM 280

AQM 350

1625 x 580 x 584 mm

1900 x 580 x 584 mm

280L

350L

2°C - 8°C

2°C - 8°C

6

6

24 Months

24 Months

10amp

10amp

Feet at front, rollers at back

Feet at front, rollers at back



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



MAINTENANCE PROGRAM

When the efficient and reliable operation of your medical refrigeration equipment is a priority, it is important to have it properly maintained.

Regular maintenance helps ensure your medical refrigerator is functioning at it's optimal performance, which in turn will extend the refrigerator's life expectancy, incur less wear and use less energy.

Our Programmed Maintenance (PM) schedule is a robust and thorough process. It is only required to be done once annually to ensure the optimum performance of your fridge.

Each PM includes:

- Full Maintenance Program
- Annual Calibration with verification certificate
- Annual Test & Tag of equipment

Our Maintenance Program is supported by the largest medical refrigeration service network in the country. Avem Quirks employs over 100 qualified Service Technicians, supported by a contractor network of 300+. We have service capability in every state and territory in Australia.

	1800 658 120
	www.avemquirks.com.au
	info@avemquirks.com.au



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

