


AQ Medical Vaccine Fridges are a **robust, cost effective** solution for all of your **vaccine refrigeration** needs.



Features

- Self closing door to ensure temp never breaches vaccine temps
- Temperature keeps from 2°C - 8°C accurately and evenly
- High-low temp audible and visual alarm
- Suitable for small vaccine bottles
- Factory installed lock ensures security
- Electronic temperature controller and display
- 24 months warranty with national service support



AQM 350

Dimensions (HxWxD)	1900 x 580 x 584 mm
Capacity	350L
Temp Range	2°C - 8°C
No. of Shelves	6
Warranty (parts and labour)	24 Months
Power Supply	10amp
Feet Included	Feet at front, rollers at back



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
 Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

