

Acapella Choice

Vibratory Positive Expiratory Pressure (PEP) System

Open airways and mobilise secretions in any patient with Acapella Choice

Numeric expiratory resistance/frequency dial
Allows therapy to be customised to meet the patients needs

One-Way Inspiration Valve
(With 22mm male fitting)
Allows inhalation without removing from the patients mouth

Patient End
Can be used with Mask or Mouthpiece



Easy to Use

- Combines the benefit of both PEP therapy and airway vibration to mobilise pulmonary secretions
- Can be used in virtually any position, patient is free to sit, stand or lay down
- Numeric dial for adjustable frequency and flow for reproducible settings
- Easy to self administer. Patient can inhale through Acapella choice without removing from mouth
- Complimentary use instruction video available upon request
- Can be used with nebulizer to provide simultaneous aerosol drug delivery while performing secretion clearance and atelectasis



Easy to Clean

Acapella Choice disassembles into four easy to clean parts. Each part can withstand autoclaving, boiling and cold sterilization. With such an easy to clean device patients can benefit from vibratory PEP therapy at home as well as in hospital.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

