

# NEW ASKIR 30

**CA-MI**

**Surgical Suction Unit**



## TECHNICAL SPECIFICATIONS

Power Feeding other voltages available upon request	230V-50hz
Max Vacuum - adjustable	-0.80 Bar / -80kPa / -600 mmHg
Max.free air flow	40 l/min
Power Consumption	107 VA
Fuse	1 x 1.6 A 250V
Intermittent duty cycle at 35°C and 110% operating voltage	120 ON / 60 OFF
Weight	Lbs 7.91 - Kg 3.60
Size	in 13.78" x 8.27" x 7.09" - mm 350 x 210 x 180
Years Warranty	2
Shipping carton	4
Manufacturing	100% in Italy

## CONFIGURATIONS

with AUTOCLAVE JAR\* with overflow valve system, antibacterial and hydrophobic filter, complete set of tubes 6 x 10mm, conical connector for probes.

\*with 1000ml autoclavable jar: RE310100/02

\*with 2000ml autoclavable jar: RE310100/03



Established 1922

**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

