

# Assut Sutures

## Polypropylene

- Ties like silk with the strength of nylon.
- Effortless passage through tissue.
- Extremely well tolerated.
- Avoids capillarity.

## AssuCryl® MonoRapid

- Absorbable monofilament synthetic made from PGCL.
- Retains 50% of strength after one week and dissolves completely within three to four months.

## AssuCryl®

- Braided, coated synthetic absorbable PGA. Predictable absorption time. High initial tensile strength.
- Easy and secure knotting. Excellent knot sliding, smooth to handle.

Available in undyed (white) and dyed (violet).



- Absorption essentially complete between 60 and 90 days.

## AssuCryl® MonoSlow Violet

- Absorbable monofilament synthetic made from PDO (polydioxanone).
- High tensile strength.
- Retains 75% of strength after three weeks and dissolves completely within six months.

## Nylon

- Blue monofilament polyamide with high tensile strength.
- The uniform and smooth surface permits effortless passage through tissue with minimum trauma.
- Extremely well tolerated.



Established 1922

## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsur (1300 633 787) Email medical@medsur.com.au Fax 07 3859 2995



# Assut sutures

Made in Switzerland

## Nylon

Blue monofilament (polyamide) non-absorbable nylon sutures.

### 45cm Thread Length <sup>3/8</sup> Reverse Cutting Needle

Order Code	Suture Gauge	Needle Size	No. in Box
AS-A5565F	2/0	24mm	36
AS-A508F	3/0	19mm	36
AS-A510F	3/0	24mm	36
AS-A505F	4/0	17mm	36
AS-A507F	4/0	19mm	36
AS-A509F	4/0	24mm	36
AS-A05-02	5/0	12mm	36
AS-A504F	5/0	17mm	36
AS-A506F	5/0	19mm	36
AS-A05-01	6/0	12mm	36
AS-A5531F	6/0	15mm	36

### 75cm Thread Length <sup>3/8</sup> Reverse Cutting Needle

Order Code	Suture Gauge	Needle Size	No. in Box
AS-A55-65	2/0	24mm	36
AS-A05-08	3/0	19mm	36
AS-A05-10	3/0	24mm	36
AS-A05-05	4/0	17mm	36
AS-A05-07	4/0	19mm	36
AS-A05-09	4/0	24mm	36
AS-A0502	5/0	12mm	36
AS-A05-04	5/0	17mm	36
AS-A05-06	5/0	19mm	36
AS-A55-31	6/0	15mm	36

## Polypropylene

Blue Monofilament with the tensile strength of a nylon suture coupled with the ease of knotting of a silk suture.

### 45cm Thread Length <sup>3/8</sup> Reverse Cutting Needle

Order Code	Suture Gauge	Needle Size	No. in Box
AS-A9565F	2/0	24mm	36
AS-A9554F	3/0	19mm	36
AS-A9564F	3/0	24mm	36
AS-A9543F	4/0	17mm	36
AS-A953F	4/0	19mm	36
AS-A9563F	4/0	24mm	36
AS-A9542F	5/0	17mm	36
AS-A9552F	5/0	19mm	36
AS-A9521F	6/0	12mm	36

### 75cm Thread Length <sup>3/8</sup> Reverse Cutting Needle

Order Code	Suture Gauge	Needle Size	No. in Box
AS-A9554	3/0	19mm	36
AS-A9564	3/0	24mm	36
AS-A9553	4/0	19mm	36
AS-A9522	5/0	12mm	36
AS-A9542	5/0	17mm	36

## AssuCryl®

Braided, coated synthetic absorbable PGA. Predictable absorption time (60 to 90 days), 50% loss after 18 days, high tensile strength with easy and secure knotting.

### 45cm Thread Length <sup>3/8</sup> Reverse Cutting Needle

Order Code	Suture Gauge	Needle Size	No. in Box	Colour
AS-A4564WF	3/0	24mm	36	Undyed
AS-A4553WF	4/0	19mm	36	Undyed

### 75cm Thread Length <sup>3/8</sup> Reverse Cutting Needle

Order Code	Suture Gauge	Needle Size	No. in Box	Colour
AS-A4544W	3/0	17mm	36	Undyed
AS-A4554	3/0	19mm	36	Violet
AS-A4564W	3/0	24mm	36	Undyed
AS-A4553	4/0	19mm	36	Violet
AS-A4563W	4/0	24mm	36	Undyed
AS-A4542W	5/0	17mm	36	Undyed
AS-A4522W	5/0	12mm	36	Undyed

## AssuCryl® MonoSlow

Monofilament synthetic, absorbable PDO (polydioxanone). Features high tensile strength, is easy to handle and passes smoothly through tissue. 25% loss after 21 days. Essentially complete within six months.

### 45cm Thread Length <sup>3/8</sup> Reverse Cutting Needle

Order Code	Suture Gauge	Needle Size	No. in Box	Colour
AS-MS4522F	5/0	12mm	36	Violet
AS-MS4521F	6/0	12mm	36	Violet

### 75cm Thread Length <sup>3/8</sup> Reverse Cutting Needle

Order Code	Suture Gauge	Needle Size	No. in Box	Colour
AS-MS4565	2/0	24mm	36	Violet
AS-MS4554	3/0	19mm	36	Violet
AS-MS4564	3/0	24mm	36	Violet
AS-MS4553	4/0	19mm	36	Violet

## AssuCryl® MonoRapid

Monofilament synthetic Polyglycolide-poly(e-caprolactone) copolymer (PGCL). Tensile strength loss of 50% at 7 days, 0% at 21 days, absorption between 90 - 120 days. Undyed.

### 45cm Thread Length <sup>3/8</sup> Reverse Cutting Needle

Order Code	Suture Gauge	Needle Size	No. in Box
AS-BMR4565WF	2/0	24mm	12
AS-BMR4554WF	3/0	19mm	12
AS-BMR4564WF	3/0	24mm	12
AS-BMR4543WF	4/0	17mm	12
AS-BMR4553WF	4/0	19mm	12
AS-BMR4522WF	5/0	12mm	12
AS-BMR4532WF	5/0	15mm	12
AS-BMR4542WF	5/0	17mm	12
AS-BMR4552WF	5/0	19mm	12



### Quality Lasts

Medical and Surgical Requisites have always supplied the most advanced technology in medicine.

This excerpt from our 1938 catalogue shows the latest innovation—an eyeless needle with thread already attached!