

# Bovie®

## SPECIALIST|PRO

*Electrosurgical Generator*



The Bovie® Specialist Pro is the electrosurgery solution for in demand physician offices and specialty facilities. Designed to operate like an operating room generator in a compact lower watt package. The multipurpose unit is ideal for monopolar and bipolar office surgical procedures.

The Specialist|PRO offers 5 common energy outputs controlled by the new sleek, push-button design, with up to 120 watts of power. The modalities include Cut, Blend, Coagulation, Fulguration and Standard Bipolar.

***Efficiency, Versatility, and Safety ...***  
*Designed benefits for you and your patients*

Industry-standard connections allow for use of a variety of monopolar and bipolar accessories providing energy versatility across a variety of specialties such as OBGYN, Urology, Plastic Surgery, Dermatology, ER and more.

Bovie NEM™ neutral electrode monitoring increases patient safety by ensuring proper pad connectivity during procedures.

The Bovie Specialist|PRO is backed by the unmatched **4-year manufacturer warranty**, offered only by Bovie.



Established 1922

**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



# BOVIE® SPECIALIST|PRO ELECTROSURGICAL GENERATOR



## STANDARD EQUIPMENT

**A1250S** electrosurgical generator,  
**ESP1** sterile disposable pencil  
**ES01** blade, **ES02** needle, **ES20** ball (1 each)  
sterile disposable electrodes  
**ESRE** (5) split disposable return electrode  
**A1252C** (2) reusable grounding cords  
110 VAC Hospital-grade power cord\*  
220 VAC Hospital-grade power cord\*  
User's guide (MC-55-237-001)\*

\* Not Pictured



**ESP1**



**ES01, ES02, ES20**



**ESRE**



**A1252C**

## COMPLEMENTARY PRODUCTS †



**SE02**  
Smoke Evacuator †



**BV-1253B**  
Monopolar Footswitch †



**A1255A**  
Footswitch Adapter †




Disposable **LLETZ/LEEP**  
Loop and Ball Electrodes †



Split Disposable Return  
Electrode with cord †



Disposable **SuperCut™**  
Needles Electrodes †



Disposable **Resistick II™**  
Electrodes †



Disposable **Dermal Loops** †

† Sold Separately

## PRODUCT SPECIFICATIONS

### DIMENSIONS Height:

15.2cm  
Width: 26cm  
Depth: 30.5cm  
Weight: < 4Kg

### OUTPUT POWER

Cut: 120 watts  
Blend: 90 watts (60% Cut, 40% Coag)  
Coagulation: 80 watts  
Fulguration: 40 watts  
Bipolar: 30 watts

### OUTPUT FREQUENCY

350 – 600 kHz

### LINE VOLTAGE

A1250S: 100-240 VAC

### LINE FREQUENCY

50 – 60 Hz

### WARRANTY

Four Years



Established 1922

## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

