



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

AMEDTEC
Medizintechnik Aue GmbH

HEALTH IS THE HEART OF OUR BUSINESS.

HOLTER-RR 24H ABPM

GENTLE AMBULATORY BLOOD PRESSURE MONITORING DIAGNOSTICS

Meaningful data: The Holter-RR guarantees reliable ambulatory blood pressure monitoring.

The Auto Feedback Logic (AFL) allows quiet and fast readings with low cuff pressure.

For precise interpretation threshold tables for children and adolescents are integrated.



MANAGE COMPLEX DATA EASILY

AMEDTEC ECGpro®

TECHNICAL DATA

Pressure measurement range	systolic 60 to 290 mmHg diastolic 30 to 195 mmHg
Accuracy	+/-3 mmHg in the range indicated
Static pressure range	0 to 300 mmHg
Pulse range	30 to 240 beats per minute
Method	oscillometric
Measurement intervals	0, 1, 2, 3, 4, 5, 6, 10, 12, 15, 20 or 30 measurements per hour
Monitoring protocols	2 programmable and 7 fixed protocols
Storage capacity	300 measurements
Operating temperatures	+10 °C bis +40 °C
Operating humidity	15 % to 90 %
Dimensions	128 mm x 75 mm x 30 mm
Weight	approx. 240 g incl. rechargeable batteries
Power supply	2 Mignon AA alkaline or rechargeable batteries
Interfaces	Cable: serial interface RS232, converter for USB Bluetooth class 1

ACCESSORIES




pouch for Holter-RR



cuffs in size XS, S, M, L, XL

AMEDTEC Medizintechnik Aue GmbH
Schneeberger Straße 5 · 08280 Aue · Germany
Phone: +49 (0)3771 59 82 70
Fax: +49 (0)3771 59 82 790
info@amedtec.de · www.amedtec.de



Manufacturer:
IEM GmbH
Cockerillstr. 69
52222 Stolberg



Certificate:
DIN EN ISO 13485

Errors and omissions excepted
M0672 018-0402-2_03