



# AQ Medical 280

## VACCINE REFRIGERATORS

### KEY FEATURES

- QCPP cold chain accredited
- Digital temperature control system preconfigured to hold temperature from +2°C to +8°C in line with 'Strive for 5' storage guidelines
- Lockable door
- Adjustable shelves
- Low energy consumption, low noise output
- Audible and visual auto resetting alarms
- Dedicated technical support line

### OPTIONAL EXTRAS

- UPS & Battery Backup
- Controller Backup Battery (temperature logging)
- AQ Box cloud temperature monitoring system
- Comprehensive programmed maintenance

### PRODUCT SPECIFICATIONS

Dimensions	Height 1575 mm
	Width 584 mm
	Depth 584 mm
Gross Weight	79 kg
Gross Capacity	280 L
Temperature Range	+2°C to +8°C
Number of shelves	6
Power Consumption/24 hrs	2.6 kWh
Refrigerant	R600a
Warranty (parts & labour)	24 months