


Integrated E-Z Split Cup II Urine Test Drug Screen

Integrated E-Z Split Cup II Urine Test Drug Screens detects DoAs in urine to the following standards:

6 Panel

- Amphetamine (300ng/ml)
- Benzodiazepine (200ng/ml)
- Cocaine (150ng/ml)
- Methamphetamine (300ng/ml)
- Morphine (300ng/ml)
- THC (50ng/ml)



Features Benefits

Six separate adulteration strips which includes oxidants/PCC, specific gravity, pH nitrite, glutaraldehyde and creatinine.

Independently validated to meet Australian Standards 4308:2023

Hinged cup lid.

Quick results.

Very strong bands.

Large flat test panel.

Urine cup is test well.

Cup lid with split key.

Key activation.

Temperature strip.

Proven history in Australia.
Extensive technical support.

Prevention of adulteration.

Lid stays attached.

High throughout.

Easy interpretation of results.

Results are easily scanned or photocopied.

Less steps for sampling.

Reduces urine exposure. Cleaner,
Sample can be saved for labs.
Prevention of adulteration.

Reliability and confidence.
Full company expertise at your disposal.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

