



Medisafe PLUS



Medisafe Plus 1000 VACCINE REFRIGERATORS

The Medisafe Plus 1000 is one of the highest quality and feature packed models of vaccine storage units available in Australia. Designed in Australia, it is purpose-built for vaccine storage in line with the Australian Government Department of Health's 'Strive for 5' vaccine storage guidelines. The Medisafe Plus 1000 is QCPP accredited giving you confidence your vaccines will be stored securely.

The Medisafe Plus range of Vaccine Refrigerators are uniquely equipped with two temperature probes for rigorous monitoring of both cabinet temperature and a simulated vaccine temperature (when using the AQ Box cloud temperature monitoring system). This yields the most accurate measurements possible and is a unique feature in the market for vaccine protection.

The AQ Box comes standard with the Medisafe Plus 1000 and includes 12 months free access to the AQ portal. The AQ box acts as an automated data logger which allows you to view and print temperature charts with Minimum and Maximum temperatures, as well as real time data from any web enabled smart device.

All Medisafe Plus units are backed by Avem Quirks' extensive service network. With qualified technicians and warehousing facilities all across Australia, you can be sure support is never far away.



Medisafe Plus 1000

VACCINE REFRIGERATORS

KEY FEATURES

- QCPP cold chain accredited
- Low energy consumption, low noise output
- Dual temperature probes for both cabinet temperature and simulated vaccine temperature
- Integrated door lock(s)
- Adjustable shelves
- Fully integrated AQ Box provides IoT cloud-based monitoring system with customizable email and text alarm notifications
- Audible and visual alarms for: temperature, door open, power failure
- Dedicated technical support line

OPTIONAL EXTRAS

- UPS Battery Backup
- Controller Backup Battery (temperature logging)
- Customisable controller parameters
- Laser printed branding (bulk orders)
- Comprehensive programmed maintenance

PRODUCT SPECIFICATIONS

Dimensions	Height: 2200mm Width: 1130mm Depth: 745mm
Gross Weight:	230 kg
Gross Capacity:	980L
Doors:	2
Temperature Range:	+2°C to +8°C
Number of shelves:	10
Power Consumption/24 hrs:	3.79kWh
Refrigerant:	R290
Warranty (parts & labour):	24 months

