

Medisafe Plus 70

VACCINE REFRIGERATORS



KEY FEATURES

- QCPP cold chain accredited
- Dual temperature probes for both cabinet temperature and simulated vaccine temperature at +2°C to +8°C
- Integrated door lock(s)
- Fully integrated AQ Box provides IoT cloud-based monitoring system with customizable email and text alarm notifications
- Audible and visual alarms for: temperature, door open, power failure

OPTIONAL EXTRAS

- UPS Battery Backup
- Controller Backup Battery (temperature logging)
- Customisable controller parameters
- Laser printed branding (bulk orders)
- Comprehensive programmed maintenance

PRODUCT SPECIFICATIONS

Dimensions	Height: 750 mm Width: 430 mm Depth: 500 mm
Net Weight:	26 kg
Gross Capacity:	70 L
Doors:	1
Temperature Range:	+2°C to +8°C
Number of shelves:	3
Power Consumption/24 hrs:	1.05 kWh
Refrigerant:	R600a
Warranty (parts & labour):	60 months



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995