



Medisafe PLUS



Medisafe Plus 420

VACCINE REFRIGERATORS

The Medisafe Plus 420 is one of the highest quality and feature packed models of vaccine storage units available in Australia. Designed in Australia, it is purpose-built for vaccine storage in line with the Australian Government Department of Health's 'Strive for 5' vaccine storage guidelines. The Medisafe Plus 420 is QCPP accredited giving you confidence your vaccines will be stored securely.

The Medisafe Plus range of Vaccine Refrigerators are uniquely equipped with two temperature probes for rigorous monitoring of both cabinet temperature and a simulated vaccine temperature (when using the AQ Box cloud temperature monitoring system). This yields the most accurate measurements possible and is a unique feature in the market for vaccine protection.

The AQ Box comes standard with the Medisafe Plus 420 and includes 12 months free access to the AQ portal. The AQ box acts as an automated data logger which allows you to view and print temperature charts with Minimum and Maximum temperatures, as well as real time data from any web enabled smart device.

All Medisafe Plus units are backed by Avem Quirks' extensive service network. With qualified technicians and warehousing facilities all across Australia, you can be sure support is never far away.



Medisafe Plus 420

VACCINE REFRIGERATORS

KEY FEATURES

- QCPP cold chain accredited
- Cutting edge new controller with Bluetooth capability, holding over a week of temperature and alarm logging
- Dual temperature probes for both cabinet temperature and simulated vaccine temperature
- Integrated door lock(s)
- 6 full size shelves and 1 half shelf
- Fully integrated AQ Box provides IoT cloud-based monitoring system with customizable email and text alarm notifications
- Audible and visual alarms for: temperature, door open, power failure

OPTIONAL EXTRAS

- UPS Battery Backup – with over 13 hours backup
- Controller Backup Battery (temperature logging)
- Customisable controller parameters
- Laser printed branding (bulk orders)
- Comprehensive programmed maintenance

PRODUCT SPECIFICATIONS

| | |
|----------------------------|--|
| Dimensions | Height: 2200mm Width: 600 mm Depth: 700 mm |
| Gross Weight: | 92 kg |
| Gross Capacity: | 422 L |
| Doors: | 1 |
| Temperature Range: | +2°C to +8°C |
| Number of shelves: | 6 |
| Power Consumption/24 hrs: | 2.46Wh |
| Refrigerant: | R290 |
| Warranty (parts & labour): | 24 months |

