

See It. Share It. Save It.



Welch Allyn Digital MacroView™ Otoscope

The Welch Allyn Digital MacroView Otoscope is a simple, affordable tool that lets you see, explain, and save digital images of the ear.

Offering you clear live and still images of the ear canal and tympanic membrane, the Digital MacroView Otoscope is a valuable tool in support of the treatment and management of your patients' ears. Easy to learn and use, the Digital MacroView is a plug and play system requiring no additional purchase of software.

Here's how the Welch Allyn Digital MacroView Otoscope can benefit you:

- > Improves communication, which can lead to better patient compliance
- > Allows you to store digital images in patient records for enhanced documentation of exam findings
- > Provides you the ability to share images with other providers as visual descriptive reference
- > Gives you an excellent training aid for educational purposes
- > Offers you the ability to view small details of tubes and hearing aids



Handle sold separately

WelchAllyn®

Advancing Frontline Care™



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Q 4077 PO Box 3965 Mt. Ommaney QLD 4074
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au www.medsurg.com.au Fax 07 3859 2995



The Welch Allyn Digital MacroView™ Otoscope

Features you expect from an otoscope and more:

USB cable for connection to your computer

Insufflation port for pneumatic otoscopy

Adjustable focusing wheel for fine-tuned views

Simple button functions
> Capture image
> Flip image
> Zoom

Digital imager for megapixel still images

Compatible with the range of Welch Allyn 3.5 V power sources

Handle sold separately

USING THE DIGITAL MACROVIEW OTOSCOPE IS EASY AS 1 – 2 – 3

1 Plug In.

Plug the otoscope into your computer using the USB cable.

2 Power On.

Put on ear specula and turn on light.

3 Play.

View live picture or capture image.

Ordering Information

23920 Digital MacroView Otoscope with 10 ft USB cable and set of four reusable ear specula

Accessories

23804 Insufflator Bulb

24400-U Additional set of four Reusable Specula

24302-U 2.5 mm Universal Poly Reusable Speculum

24303-U 3.0 mm Universal Poly Reusable Speculum

24304-U 4.0 mm Universal Poly Reusable Speculum

24305-U 5.0 mm Universal Poly Reusable Speculum

24320 SofSeal Small, Disposable Sealing Tip

24330 SofSeal Medium, Disposable Sealing Tip

52432-U 2.75 mm Single-use Ear Specula

52434-U 4.25 mm Single-use Ear Specula

06500-U Lamp

23900-10 Replacement USB Cable, 10 ft

71000-A NiCad Rechargeable Handle (Plug-in)

71000-C NiCad/Convertible Rechargeable Handle (Plug-in)

71670 NiCad Rechargeable Handle for Universal Charger

71900 Lithium Ion Rechargeable Handle with AC Charging Module

71910 Lithium Ion Rechargeable Handle for Universal Charger

System Requirements: Windows® XP Service Package 2 or higher with an available USB 2.0 port.

To find out how the new Welch Allyn Digital MacroView Otoscope can help you see, explain and save digital images of the ear, go to www.welchallyn.com or call 1.800.535.6663.

4341 State Street Road, Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361 www.welchallyn.com

WelchAllyn®

Advancing Frontline Care™

© 2008 Welch Allyn SM2978 Rev C



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Q 4077 PO Box 3965 Mt. Ommaney QLD 4074
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au www.medsurg.com.au Fax 07 3859 2995



Established 1922