

CRYOSPRAY

PROBE REPLACEMENT GUIDE

Cryospray
Model No.
WH-099



Models manufactured prior to **30th November 2011**
with Serial Numbers up to **12111279**

Use Cryospray Probe Part No: **31-CP099**



31-CP099

Cryospray
Model No.
WH-121



Cryospray Units manufactured from **1st December 2012** to **30th May 2016**
Serial Numbers: **12112180 – 1604111-25**

Use Cryospray Probe Part No: **31-CP121**



31-CP121

Serial No. Locations:

- A) Bottom of Flask
- B) Under Turret

A



B

