

LIQUID NITROGEN DEWAR'S

PRICELIST - MARCH 2017

Designed specifically for liquid nitrogen storage, Walmay Liquid Nitrogen Dewar's are characterised by their small calibers and low static evaporation loss.

- High strength aluminium alloy structure
- Multi layered thermal isolation design
- Protective jacket to prevent collision

WMPYDS-10L	1-5 Units
	6+ Units
WMPYDS-20L	1-5 Units
	6+ Units

Above prices are per unit and exclude GST and delivery.

Free level
measuring
scale with every
Walmay Dewar



Model	WMPYDS-10L	WMPYDS-20L
Part Number*	20601714	20601711
Capacity	10.0 L	20.0 L
Empty Weight	6.2 kg	11.2 kg
Neck Diameter	50 mm	50 mm
Outer Diameter	303 mm	409 mm
Height	530 mm	650 mm
Static Evaporation Loss**	0.10 L/d	0.10 L/d
Static Hold Time	101 d	202 d

* Containers with same PN may contain various configurations. Please confirm when ordering.

** Static evaporation loss and static hold time are nominal. Actual rate and hold time will be affected by the nature of container use, atmospheric conditions, and manufacturing tolerances.