



CAM-RE300200 is a piston compressor nebulizer for aerosoltherapy. Indicated for home-care and intensive use, its design includes a convenient handle for easy transport of unit and a large compartment for storage of accessories and medications.

INCLUDED ACCESSORIES

HI-FLO Jet Nebulizer
Adult Masks with elastic strings
Child Masks with elastic strings
Mouthpiece
Nosepiece (non invasive)
Air Tube
Air Filter

FEATURES

Compressor	Oiless and maintenance-free piston pump
Power Feeding	230 V-50 Hz
Power consumption	170 VA
Max Pressure	3,0 bar
Operating Pressure	1,3 bar 18,9 psi 130 Kpa
Max Air Flow	16 L/min
Operating Air Flow	5,2 L/min
Neb-rate with 4 ml solution	0,5 ml/min
Mass Median Aerodynamic Diameter	3,25 μm (MMAD)
Respirable Fraction (< 5 μm)	> 70 %
Noise Level (EN 13544-1 Standard)	60 dB
Weight	1,60 Kg
Size (Main Unit)	cm 13,8x11,7x23,5
Years of Warranty	4
Place of Manufacturing	Italy



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

