

CA-MI Flo Eolo Nebuliser



The Eolo is a piston compressor nebulizer (nebuliser) which is for aerosol therapy. Designed for home care with convenient handle and a large compartment for easy transport of unit, accessories and medication. This product features recently upgraded for non-stop operation in case intensive treatments are required.

Included with the nebuliser -

- HI-FLO Jet Nebulizer
- Adult Nebuliser Mask
- Pediatric Nebuliser Mask
- Nebuliser Mouthpiece
- Air Tube
- Air Filter

Importnat Nebulizer features -

- Max Pressure 2.5 Bar 36 psi 250 kPa
- Operating Pressure 1.1 Bar 16 psi 110 kPa
- Max Air Flow 5 l/min
- Neb-rate with 4 ml solution 0.40 ml/ min
- Mass Median Aerodynamic Diameter (MMAD) according to EN- 13544-1 2,44 ?m
- Respirable Fraction (< 5 ?m) > 80%
- Duty Cycle Non- stop operation
- Noise Level (EN 13544-1 Standard) Approx. 55 dBA
- Weight 1,65 Kg



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

