


DBD-167-19A Clonal

Drug Screen Test Cup

Integrated E-Z Split Cup II Urine Test Drug Screens detect DoAs in urine to the following standards:

Amphetamine (300 ng/ml),
Benzodiazepine (200 ng/ml),
Cocaine (300 ng/ml),
Methamphetamine (300 ng/ml),
Morphine (300 ng/ml),
and THC (50 ng/ml).



Features

Six separate adulteration strips which includes oxidants/PCC, specific gravity, pH, nitrite, glutaraldehyde and creatinine.

Hinged cup lid

Independently validated to meet Australian Standards 4308.

Quick results.

Very strong bands.

Large flat Test Panel.

Urine cup is test well.

Cup lid with split key.

Key activation.

Temperature strip.

Proven history in Australia.

Extensive technical support.

Benefits

Prevention of adulteration.

Lid stays attached.

Consistent with standards for cut off levels.

High throughput.

Easy interpretation of results.

Results are easily scanned or photocopied.

Less steps for sampling.

Reduces urine exposure, cleaner. Sample can be saved for labs.

Prevention of adulteration.

Reliability and confidence.

Full company expertise at your disposal.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

