

# Introducing the CoaguChek® Pro II

*The next generation that makes all the difference  
for compliance and effective coagulation monitoring*



When it comes to patient management, the right tool can make all the difference. New CoaguChek Pro II brings simplicity, accuracy, and connectivity together for coagulation monitoring that will impact both you and your patients.

► Simple ► Fast ► Easy to Use ► Reliable

**CoaguChek®**



Established 1922

**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



# Designed for immediate use

The new CoaguChek Pro II allows you to make a seamless transition to the next generation in patient anticoagulation monitoring.



## Ease of use and convenience

- Easy-to-use intuitive interface
- Touch screen operation
- Minimum 8µL drop of capillary, venous or arterial fresh blood
- Out of meter side and top dosing on the strip
- Results within minutes for PT/INR and aPTT
- Large memory capacity with 2000 patient and 500 QC results



## Quality and reliable results

- On-board single-channel strip control (OS2C) checking to automatically ensure the strip's integrity
- Precise results, correlating very well to laboratory methods<sup>1,2</sup>
- Increase patient's compliance and time in therapeutic range (TTR)<sup>3</sup>
- Insensitive to heparin allowing for the rapid transition from heparin to warfarin therapy
- Better patient outcomes with lower morbidity<sup>3</sup>



## Workflow efficiency for CoaguChek Pro II

- With built-in WiFi to enable WLAN connectivity
- 2-D barcode scanner
- QR code communication capability
- Patient/operator ID management and printer option
- Quality control (QC) and operator lockouts

## Creates highly efficient workflow

The new CoaguChek Pro II puts advanced PoC coagulation testing right into your practice.



### Barcode scanning

Frees up manual time for clinicians and patients



### Accurate results in 1 minute<sup>4</sup>

Confidence from the start



### Automatic wireless data transfer

Ensures data integrity and immediately accessible results



### Clinicians can take action sooner

Optimises treatment efficiency



### Supports better outcomes<sup>4,5,6</sup>

From patient health to overall efficiency

New CoaguChek Pro II makes compliance easy, every step of the way



Established 1922


## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



# Introducing CoaguChek Pro II

*Realise significant benefits at every touchpoint*



**Contact your Roche representative**  
and put it to work for you and your patients



Established 1922

**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
Phone 1300medsurg (1300 633 787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995



**More than 2 million patients on vitamin K antagonists around the world are testing with CoaguChek® systems,\* so you can feel confident that you are making the right decision**

**References:**

1. Plesch W et al. Results of the performance verification of the CoaguChek XS system. *Thromb Res* 2008. 123: 381-389.
  2. Plesch WJ et al. CoaguChek XS Plus system are equivalent to those of the XS system. *J Thromb Haemost* 2007. 5 Supplement 2: P-T-108.
  3. Chiquette E., Amato M., Bussey H. "Comparison of an Anticoagulation Clinic With Usual Medical Care." *Arch Intern Med*. 1998; 158:1641-1647.
  4. CoaguChek® PT Test [Method Sheet]. Rotkreuz, Switzerland: Roche; 2015.
  5. Ward A, Tompson A, Fitzmaurice D, Sutton S, Perera R, Heneghan C. Cohort study of Anticoagulation Self-Monitoring (CASM): a prospective study of its effectiveness in the community. *Br J Gen Pract*. 2015;65(636):e428-437.
  6. Khan TI, Kamali F, Kesteven P, Avery P, Wynne H. The value of education and self-monitoring in the management of warfarin therapy in older patients with unstable control of anticoagulation. *Br J Haematol*. 2004;126(4):557-654
- \* Data on file

COAGUCHEK is a trademark of Roche.

© 2016 Roche.

Roche Diagnostics Australia Pty Limited  
 ABN 29 003 001 205  
 PO Box 955  
 Castle Hill NSW 1765 Australia  
 Phone: 02 9860 2222  
[www.coaguheck.com](http://www.coaguheck.com)

**CoaguChek®**



**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
 Phone 1300medsurg (1300 633 787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995

