

SP-20

Handheld Pulse Oximeter



Features

- Measuring SpO2, Pulse Rate and Temperature(optional)
- 3.5 inch TFT display of PI%, plethysmograph and pulse bar
- Adult, pediatric and neonate measurement
- Lithium battery up to 18 hours
- 500 hours SpO2 data storage
- PC software for data upload

Technical Specifications

SpO2

Technique: Dual-wavelength LED
Measuring range: 35%~100%
Measuring accuracy: $\pm 3\%$ (range from 70% to 100%)

Pulse Rate

Measuring range: 30bpm~240bpm
Measuring accuracy: ± 2 bpm or $\pm 2\%$

Temperature(Optional)

Measuring range: 32 C ~ 43 C
Measuring accuracy: ± 0.2 C
Measuring time: < 5 seconds
Measuring probe: Infrared ear probe

Measuring Mode

Continuous or Spot check mode

Perfusion Index

PI Range: 0.2%~20%

Data Storage

Up to 500 hours

Display

3.5 TFT display

Power Supply

Lithium battery

Charging

Charging stand or USB cable

Dimension

158 mm x 73 mm x 25 mm

Net Weight

230g

Data Upload

USB cable or wireless

Operating Environment

Operating Temperature: 5°C ~40°C
Operating Humidity: 15%-93%



Adult clip sensor



Adult soft sensor



Pediatric clip sensor



Pediatric soft sensor



Neonate wrap sensor



Infrared ear thermometer



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

