

Dalcross Surgeon's Stools



Surgeon Stool with Footrest

Code: 1381

Surgeon Stool with Backrest & Footrest

Code: 1383

- Height adjustment 610mm to 870mm
- Hand control lever
- 90mm Thick high-density foam padded seat
- Seat diameter 380mm
- Vinyl upholstery colours: black or grey
- Safe working load: 140kg

Surgeon Stool

Code: 1380

Surgeon Stool with Backrest

Code: 1382

- Height adjustment 530mm to 720mm
- Hand control lever
- 90mm Thick high-density foam padded seat
- Seat diameter 380mm
- Vinyl upholstery colours: black or grey
- Safe working load: 140kg

Poly Top Surgeon Stool

Code: 1380P

- Height adjustment 470mm to 660mm
- Hand control lever
- 35cm Diameter soft polyurethane seat
- Durable, withstands impacts and sharp objects
- Safe working load: 140kg



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

