

Dalcross Electric Couches



1 Motor 2 Section

Code: 1242M

- Footswitch controlled electric motor for Hi-Lo lift function
- Gas strut assisted backrest
- Minimum height: 460mm
- Maximum height: 910mm
- Overall dimensions 1955mmL x 710mmW
- 200kg Safe working load
- Quality upholstery finish in grey or blue vinyl

1 Motor 3 Section

Code: 1243M

- Footswitch controlled electric motor for Hi-Lo lift function
- Gas strut assisted backrest and legrest
This means no effort to raise (without a patient)
- Minimum height: 460mm
- Maximum height: 910mm
- Overall dimensions: 1955mmL x 710mmW
- 200kg Safe working load
- Quality upholstery finish in grey or blue vinyl

All-Electric 2 Section

Code: 1242E

- Handswitch controlled electric motors for Hi-Lo lift function and backrest
- Optional dual control kit available - Footswitch for Hi-Lo lift, Handswitch for back
- Minimum height: 460mm
- Maximum height: 910mm
- Overall dimensions: 1955mmL x 710mmW
- 200kg Safe working load
- Quality upholstery finish in grey or blue vinyl

All-Electric 3 Section

Code: 1243E

- Handswitch controlled electric motors for Hi-Lo lift function, backrest and legrest
- Optional dual control kit available - Footswitch for Hi-Lo lift, Handswitch for backrest and legrest control
- Minimum height: 460mm
- Maximum height: 910mm
- Overall dimensions: 1955mmL x 710mmW
- 200kg Safe working load
- Quality upholstery finish in grey or blue vinyl



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



Dalcross Electric Couches



Bariatric 2 Section

Code: 1252

- Handswitch controlled electric motors for Hi-Lo lift function and backrest
- Optional dual control kit available - Footswitch for Hi-Lo lift, Handswitch for backrest
- Minimum height: 460mm
- Maximum height: 910mm
- Overall dimensions: 1955mmL x 780mmW
- 330kg Safe working load
- Quality upholstery finish in blue vinyl



Bariatric 3 Section

Code: 1253

- Handswitch controlled electric motors for Hi-Lo lift function, backrest and legrest
- Optional dual control kit available - Footswitch for Hi-Lo lift, Handswitch for backrest and legrest control
- Minimum height: 460mm
- Maximum height: 910mm
- Overall dimensions: 1955mmL x 780mmW
- 200kg Safe working load
- Quality upholstery finish in blue vinyl



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

