

Dalcross Steel Frame Couches



Standard Exam Couch with Powdercoated Frame

Code: 1001

- Powdercoated frame colours: grey or beige
- Vinyl upholstery colours: grey, black, blue or beige
- Flat packed for self assembly
- Adjustable headrest: 610mmL x 610mmW
- Overall dimensions: 1820mmL x 610mmW x 780mmH
- Safe working load 125kg



Heavy Duty Exam Couch

Code: 1004

- Grey powdercoated frame
- Vinyl upholstery colours: grey or black
- Comfortable thick foam padding
- Flat packed for self assembly
- Adjustable headrest: 610mmL x 610mmW
- Overall dimensions: 1820mmL x 610mmW x 780mmH
- Safe working load 250kg

Optional bolt-on metal sided drawer set:

- Model 1048 Drawers set of 3 for 1004 couch



Mobile Examination Couch

Code: 1005

- Sturdy chrome finish tubular frame
- Vinyl upholstery colours: grey or black
- Fixed 75mm thick foam padded top
- 125mm Rubber castors, 2 with locks
- Flat packed for self assembly
- Adjustable headrest: 610mm x 610mmW
- Overall dimensions: 1820mmL x 610mmW x 780mmH
- Safe working load: 125kg

Add accessories to form a basic patient trolley:

- Safety rails
- wire basket
- IV pole & clamp



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

