

Dalcross Timber Combination Couch



Timber Cabinet Couch

Code: 1203

- Cabinet made from colourboard with impact resistant pvc edging
- Sturdy construction and quality finish
- Cabinet colours: lemon tree, beech, jarrah, grey
- Vinyl upholstery colours: windsor blue, charcoal, black or doeskin
- Custom cabinet and vinyl colours available
- Five drawers
- Sliding panel doors with inner shelf
- Retractable instrument shelf
- Adjustable headrest: 580mmL x 585mmW
- Overall dimensions: 1790mmL x 585mmW x 820mmH
- Safe working load 150kg
- Matching single and double step-up stools available

Vinyl Colours



Charcoal

Black



Doeskin

Windsor

Cabinet Colours



Light Grey

Jarrah

Beech

Lemon Tree



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

