

# DermLite DL4W



## The most refined vision in dermoscopy.

The newest 4th generation DermLite achieves brighter white illumination in both the polarized and non-polarized modes. With its amazingly smooth, ergonomically optimized all-aluminum design that comfortably molds to your hand, DermLite DL4W puts phenomenal versatility, the largest optics of any pocket-size DermLite, and instant on/off control right at your fingertips.

## Larger lens, smaller package.

The DL4W has substantially brighter polarized & non-polarized white illumination when compared to the DL4 model, features a large field of view, and is more than 25% smaller than the previous generation DL3.

## Quick workflows.

DL4 uses snap connectivity to achieve groundbreaking convenience and fast workflows. Not only can the glass faceplate (with 10 mm reticle) simply snapped on and off, your mobile devices – iPad, iPhone or other smartphones – can attach magnetically, enabling you to go between taking clinical and dermoscopic photos quicker than ever before.



Established 1922

## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



## **Infection Control? You're covered.**

Along with the growing need for dermoscopy, infection control is becoming an increasingly important aspect of providing care. With its Infection Control and Epidemiology system IceCap™, DL4W has you covered. Disposable snap-on IceCaps attach to the spacer and may greatly reduce the risk of cross-contamination.

## **Smooth. Speedy. Precise.**

A newly developed DL4W retractable glass faceplate design allows you to go between non-contact & contact modalities smoother than ever and lets you precisely focus your image while in skin contact.

## **See deeper. Even into the battery.**

DermLite DL4W's high-powered LEDs and rechargeable battery produce light in an extremely energy-efficient way. And, in order to know when you run out of battery, DL4W now has a more granular four-level charge indicator. If you need to recharge, you can do so from any USB port using the standard Micro USB to USB cable included with each DL4W.

In its elegance and compactness, DermLite DL4W is simply the most capable dermatoscope available. Includes DermLite DL4W device, USB charger, Napa Leather pouch, silicone sleeve, lens cloth, and 5 IceCaps.

Warranty: 10 years



Established 1922

## **Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
Phone 1300medsurg (1300 633 787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995

