

DermLite Lumio®



Lumio is the ideal device for general skin exams where an in-depth, glare-free image of the skin is essential.

With its 75 mm lens offering a magnification of 2x, the power of 40 bright-white LEDs (light emitting diodes), and cross-polarization, Lumio is a great tool for a variety of dermatological exams: Varicose veins, pigmented skin lesions, hair follicles and much more.

Lumio is powered by AA batteries that allow for a minimum of 250 exams at 30 seconds each.

Model No. LUM is supplied with one Lumio, four AA batteries and one protective neoprene pouch.

SNAPSHOT
permanent cross-polarization
40 LEDs
AA batteries
75 mm 2x lens
stylish, ergonomic design
protective neoprene pouch

Warranty
5 years



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

