

DermLite DL100



FIRST IN NO-OIL POCKET EPILUMINESCENCE. SIMPLY INGENIOUS.

Introducing DermLite DL100, the original choice for oil-free pocket epiluminescence microscopy. DermLite DL100 speeds up the early detection of melanoma with enhanced visualization of pigmented skin lesions - without the use of oil.

By incorporating diode lighting technology for bright-white, natural illumination, a 10x Hastings Triplet lens and advanced cross-polarization technology, DermLite allows for a clear and deep view into pigmented lesions for early detection of melanoma and dysplastic nevi in a matter of seconds.

Enhanced Visualization of Basal Cell Carcinoma (BCC)

With 3Gen's signature [dry] epiluminescence and the elimination of a glass faceplate, visualization of telangiectasia associated with the early detection of BCC is greatly enhanced. Vascular patterns can be seen more easily with DermLite compared to other hand-held dermatoscopes, since the vessels are not missed because of the glass-slide compression. (Dr. Guiseppe Argenziano, Naples, Italy).



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



A Versatile Powerhouse

The versatile DermLite has a multitude of applications. Viewing pigmented lesions, varicose veins, hair transplants, skin irritations, or skin lesions is just the beginning. Together with the DermLite's compact size and lightweight design (only 3.5 oz/99g), the possibilities are virtually infinite. To make sure you will not run out of power, the DermLite's endurance measures over 1200 lesion examinations* per battery, due to its diode lighting technology and long-life lithium power.

* at 30 sec/lesion



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

