

Lumio S



Lumio®S opens a new era in skin diagnostics. With its superb 40 mm lens producing 4x magnification, this convenient device hits the sweet spot between large-area, low magnification viewers and the close quarters of a dermoscopic image.

In fact, Lumio S lets you see an area that is 4 times larger than what your dermatoscope provides, while showing lesion structures invisible to the naked eye. Simply push & hold to turn on the device, quick-tap to change between two sets of polarized and non-polarized LEDs.

The minimalist design of Lumio S is even more functional than simply cool. A clean, elegant exterior with only a single USB jack is much easier to keep clean between patients, its ergonomic soft-touch body provides a comfortable, secure grip, and its genuine Corning®Gorilla®Glass top and bottom surfaces are extremely strong and scratch-resistant.

Lumio S may be charged from any USB power source using the included USB cable.

Key features:

- Polarized & non-polarized
- Large field-of-view
- See lesion structures
- Charge via USB



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

