

Edan



Edan SE-601

- Large backlit colour screen
- Multi channel (3,4,6) dual power supply
- Fully interpretive
- 200 ECG memory

SE-601

Acquisition Mode	Simultaneous 12 lead acquisition (Defibrillation protection built-in)
Recording Mode	Auto / Manual / Rhythm
Rhythm Lead	Single lead or three leads selectable
Sensitivity	2.5, 5, 10, 20, 10/5, AGC (mm/mV)
Speed	5mm/s, 10mm/s, 12.5mm/s, 25mm/s, 50mm/s (± 3%)
Filter	EMG filter 25Hz/35Hz/45Hz/OFF Anti Baseline Drift - Anti-baseline drift 0.05/0.15/0.25/0.5/0.32/0.67Hz DFT filter 0.05/0.15/0.25/0.5/0.32/0.67Hz Low pass filter 150Hz/100Hz/75Hz AC filter 50Hz/60Hz
Acquisition technology	24 Bits A/D converter - CMRR >115dB Sampling rate 1kHz/channel, 10,000Hz for pacemaker detection
Display	5.7" 640 x 480 TFT colour LCD
Printer	Thermal printer - 110mm*140mm*144P (Z-fold)
Power Supply	AC: 100V-240V, 50/60Hz
Battery	Built in rechargeable Li-ion battery (approx. printing 280 ECG)
Communication ports	Ethernet, RS232, USB, (WIFI optional)
Safety Standards	IEC 60601-1: 1988+A1+A2, EN 60601-1:1990+A1+A2, IEC/EN60601-2: 2001+A1, IEC/EN60601-2-25, ANSI/AAMI EC-11
Certificate	Australia TGA: ARTG 178086 Class IIa
Dimension	32.2cmx31cmx10.1cm
Weight	approximately 3kg (excluding recording paper and battery)



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

