

Connecting



from

heart

to

heart



SonoTrax Series Ultrasonic Pocket Doppler can meet the routine examination requirements of obstetricians, midwives and pregnant women.

SonoTrax

Ultrasonic Pocket Doppler

Features

- Easy-to-use
- Ergonomic design
- Interchangeable probe
- Built-in audio recorder
- Large backlight LCD display
- Accurate FHR detection with clear sound
- Support 2MHz, 3MHz, 4MHz, 8MHz waterproof probes
- Automatic and manual FHR counting
- Long time continuously use (4 hours)
- Automatic power-off



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



SonoTrax

Ultrasonic Pocket Doppler



Specifications

- Display:** 45x25mm LCD display
- Heart rate range:** 50~210bpm
- Resolution:** 1bpm
- Audio sampling frequency:** 4KHz /8KHz
- Recording length:** 480s/240s
- Audio output power:** 0.5W
- Auto shut-off:** 1 minute
- Size:** 32 (Depth) x 85 (Width) x 138 (Height) mm
- Weight:** 290±5g (including one battery)

Model	Interchangeable probes	FHR display	Built-in audio recorder	Backlight	Waterproof probe	Probe auto-detection display	Battery status display	Low battery indication	AC recharge
SonoTrax Lite	Yes	No	No	No	Yes	No	No	Yes	No
SonoTrax Basic	Yes	Yes	No	No	Yes	Yes	No	Yes	No
SonoTrax Basic A	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes	No
SonoTrax Pro	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No
SonoTrax II	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes	Yes
SonoTrax II Pro	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Optional accessories	Rechargeable Li-ion battery, charger, earphone, carry bag, 9V alkaline battery, adapter								



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
 Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

