

Multigate Epidural Pack



Product code: MM-24-023

Pack Contents

Description	Size	Quantity in Pack
Sterile Field Blue	90cm x 90cm	1
Drape 150cm x 95cm, Square Aperture	Medium	1
Tray with 4 compartments Blue	23.7cm x 15.3cm x 2.3cm	1
Gallipot Graduated	20-120mL	1
Gauze Swabs Green	10cm x 10cm	5
Ramsey Forceps - Plastic Blue	23cm	1
Plastic Forceps Yellow	107mm	1
Syringe Luer Slip	5mL	1
Syringe Luer Slip Concentric	2-3mL	1
Needle 27G x 1/2" x 12mm		1
Needle 18G x 1 1/2" Blunt Draw Up		1



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

