

TREATMENT TABLES

Three Section Treatment Table

Ideal for Physiotherapy, Massage, Osteopathy and Beauty Treatments.

Every Practice should have at least one heavy duty table allowing the safe lift of heavy patients reducing fatigue and back strain on the Practitioner.

The manual adjustable head section includes a nose hole and can be lowered allowing the prone patient to fully relax the neck muscles.

The manual adjustable leg rest can be raised enabling the table to turn into a chair for when you need your patients in a seated position.



Features

- Electric height adjustment with handswitch
- Column lift eliminating total crush points
- Independently lockable castors
- Durable powder coated rust resistant frame
- High quality anti-microbial vinyl & CFC free foam
- 5 year warranty on frame & electrics
- Australian designed & manufactured

Specifications

- 300kg Maximum Lifting Capacity (moving)
- 400kg Maximum Working Load (stationary)
- 470mm Minimum Height
- 870mm Maximum Height
- 610mm Table Width
- 1840mm Table Length
- Manual Headrest adjustment +40° to -90°
- Manual Leg Rest Adjustment +70°
- Nose Hole in headrest as standard



The Healthtec Warranty

All Healthtec™ manufactured items are covered for defects in material or workmanship for 5 years from date of purchase. Original suppliers warranty applies.

TGA AR-00106892-6 © October 2013






Supporting you

www.healthtec.com.au

PRICE SHEET

Three Section Treatment Table

	CODE	QTY
Product Description		
Evolution2 Treatment Table (610 wide) – Three Section w/castors	55251	
Width Options		
670mm Width Option	-65	
710mm Width Option	-700	
770mm Width Option	-75	
Custom Width Option	6014	
Foam & Vinyl Options		
75mm Comfort Foam	6005	
Non-Stock Vinyl Option	6018	
Leather Upgrade	6008	
*please refer to our colour chart for vinyl colour options		
Nose Hole Options		
Nose Hole – Extra (Footend)	6011	
Nose Hole Plug	6010	
Face Crest Pro (Footend)	5519	
Armrest & Footrest Options		
Procedure Arm	5102	
Prone Armrest – Massage	5105-3	
Prone Armrest (for facecrest pro)	5119	
Footrest Extension	5104	
Gynae Leg Supports	5114	
Table Protection		
Paper Towel Roll Holder	6017	
Fitted Head Rest Cover (for Face Crest Pro) 100 per pack	ANWRC	
Fitted Treatment Table Cover w/nose hole (Micro Fabric)	7201	
Extended Warranty		
10 years Frame & Electrics	7800	
Freight		
Freight charges are not included in our prices. Please contact your Healthtec representative for freight costs and options.		
Shipping Dimensions	2050mm x 850mm x 650mm	
Shipping Weight	80kg (approximately)	



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

