

EXAMINATION TABLES

Examination Table

Every practice should have at least one heavy duty table allowing the safe lift of heavy patients reducing fatigue and back strain on the practitioner.

The electric HI-LO can safely lift up to 300kg, while the low minimum height allows for easy patient access and safe wheelchair transfers.

The manual adjusting backrest (electric upgrade available) enables a comfortable seated or semi reclined position, while the column profile shortens the distance between practitioner and patient to provide unmatched ergonomic accessibility while providing the greatest stability, comfort and portability.



Features

- Electric Height Adjustment with handswitch
- Column lift eliminating crush points
- Independently Lockable Castors
- Durable powder coated rust resistant frame
- High quality anti-microbial vinyl & CFC free foam
- 5 year warranty on frame & electrics
- Australian designed & manufactured

Specifications

- 300kg Maximum Lifting Capacity (moving)
- 400kg Maximum Working Load (stationary)
- 470mm Minimum Height
- 870mm Maximum Height
- 710mm Table Width
- 1930mm Table Length
- Manual backrest adjustment +70° (electric upgrade available)

Available options:



5102 – Procedure Arm



5114 – Gynae Leg Supports



5112 – Side Rails



6017 – Paper Towel Roll Holder



The Healthtec Warranty

All Healthtec™ manufactured items are covered for defects in material or workmanship for 5 years from date of purchase. Original suppliers warranty applies.

TGA AR-00106892-6

© September 2013



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



PRICE SHEET

Examination Table

Product Description	CODE	QTY
Product Description		
Evolution2 Examination Table (710 wide) w/castors	55601	
Evolution2 Examination Table (710 wide) w/Electric backrest & castors	55601-EB	
Width Options		
750mm Width Option	-75	
Custom Width Option	60140	
Foam & Vinyl Options		
75mm Comfort Foam	6005	
Non-Stock Vinyl Option	6018	
Leather Upgrade	6008	
*please refer to our colour chart for vinyl colour options		
Nose Hole Options		
Nose Hole – Examination Backrest	6011-2	
Nose Hole Plug	6010	
Examination Options		
Procedure Arm	5102	
Gynae Leg Supports	51140	
Side Rail (Left/Right) Each	5115	
Drip Pole (IV) Mounting Bracket	6041	
Drawer Options		
Side Drawer	5120	
End Drawer	5121	
Table Protection		
Paper Towel Roll Holder	6017	
Disposable Fitted Bed Sheet – 10 per pac	ANWFIT	
Extended Warranty		
10 years Frame & Electrics	7800	
Freight		
Freight charges are not included in our prices. Please contact your Healthtec representative for freight costs and options.		
Shipping Dimensions	2050mm x 850mm x 650mm	
Shipping Weight	80kg (approximately)	

