

FLEXLED 15 Exam Light

Powerful. Adjustable. Inexpensive.

Ideal for any general practitioner, specialist and hospital examination including gynae and dermatology use.



FlexLED with Wall Mount
Model 1500



FlexLED on Mobile Base
Model 1501

- Seven powerful long-life LED's with focused beam
- LED's provide 50,000 Lux @ 0.5m
- 1000mm Long flexible gooseneck arm
- Variable intensity and colour control
- High energy efficiency
- Low heat and no infrared radiation
- Low maintenance. No globes to change.
- Mounting options: wall or mobile floor stand

Touchpad Control on back of lamphead



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

