

BETA 200® Diagnostic Sets



The BETA 200 Diagnostic Set includes a BETA 200 Ophthalmoscope, BETA 200 Fibre Optic Otoscope, one set of reusable tips, ten disposable tips and one spare bulb each.

BETA 200 Fibre Optic Otoscope:

- All-metal housing matt-chrome plated
- Instrument head matt-black inside
- Distal Fibre Optic illumination with XHL Xenon Halogen Technology
- Glass viewing window made of multi-coated optical glass with 3x magnification
- Swivelling viewing window built into instrument
- Insufflation port

BETA 200 Ophthalmoscope

- 27 lenses from -35D to +40D
- Eliminates corneal reflex
- Ideal for small pupils
- Dustproof



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

