

## GAMMA Sphygmomanometers



### GAMMA G7®

- Shockproof
- Latex-free
- Exclusive pressure-distribution system
- Pressure-proof manometer
- Micro-filter
- Large-diameter scale
- High quality indicator
- Large insufflation bulb
- Exclusive air-release valve
- Upper housing made of aluminium alloy
- Lower housing is moulded of a 2 component thermo plastic (Polyamid/TP Elastomer)
- Light alloy, adjustable spoon
- High quality, single-tube cuff with velcro fastener
- Ergonomic design
- 10 year warranty\*

\*against defects in manufacture and workmanship for a period of 10 years from the date of purchase

### GAMMA G5®

Combines all the design and quality features of the GAMMA G7, except:

- Housing is moulded of a 2 component thermo plastic
- Adjustable spoon (Polyamide thermo plastic)
- 5 year warranty\*

\*against defects in manufacture and workmanship for a period of 5 years from the date of purchase

### GAMMA GP®

Combines all the design and quality features of the GAMMA G5 but with a push-button valve:

- Innovative push-button valve
- 5 year warranty\*

\*against defects in manufacture and workmanship for a period of 5 years from the date of purchase



Established 1922

## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



## GAMMA Range

### GAMMA XXL LF®

- Latex-free
- Durable polycarbonate construction
- Oversized 130mm diameter scale
- Long spiral inflation tube (3m)
- High quality indicator
- Precision air-release valve
- Pressure-proof manometer
- Wall, desk, rail and wheeled stand
- 10 year warranty\*

\*against defects in manufacture and workmanship for a period of 10 years from the date of purchase



Established 1922

## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900  
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077  
Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au)

[www.medsurg.com.au](http://www.medsurg.com.au)  
Fax 07 3859 2995

