

K180® Diagnostic Sets

The K180 Diagnostic Set includes a K180 Ophthalmoscope, K180 Fibre Optic Otoscope and one set of reusable tips.

K180 Fibre Optic Otoscope:

- Fibre Optic illumination
- Polycarbonate construction
- 3 x magnification scratch-resistant acrylic viewing window
- Otherwise technically identical to the BETA 200 Fibre Optic Otoscope

K180 Ophthalmoscope

- Precision spherical HEINE optics
- XHL Xenon Halogen Technology
- 27 lenses from -35D to +40D
- Ergonomical shape
- Soft orbital rest
- Choice of 5 apertures
- Dust-resistant housing
- High-impact polycarbonate housing for durability



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

