

## ∴ HEINE BETA® 400 F.O. Otoscope

The highest magnification for the best diagnosis.

NEW!

4.2x

MAGNIFICATION



DEPTH OF FIELD



ONE-HAND  
OPERATION

THE HEINE BETA® 400 FIBER OPTIC OTOSCOPE.  
4.2x MAGNIFICATION. EXCEPTIONAL DEPTH OF FIELD. ONE HAND OPERATION.

The HEINE BETA 400 F.O. Otoscope is the most advanced high quality otoscope making more details visible, while maintaining familiar handling and functionality. With a previously unmatched magnification factor of 4.2x and multi-coated precision optics for sharp spatial images, it truly provides for the best diagnosis. Combined with a large field of view the ear canal and tympanum are clearly visible at a glance. The unique opening mechanism (patent pending) can be operated with one hand and flips up the optical system. Application of instruments through the device becomes easily possible.



**4.2x**

**4.2x MAGNIFICATION**

Unrivalled: More detail of anatomical structures and easy recognition of foreign bodies.



**ONE HAND OPERATION**

Easy thumb control: Optical system to be flipped open and out of the way, enabling instrumentation and ease of cleaning.



**EXCEPTIONAL DEPTH OF FIELD**

Precise image: Provides sharp spatial images in a convenient working range.

BETA 400 F.O. Otoscope **without handle** [01]

**2.5V: B-001.11.400**

**3.5V: B-002.11.400**

BETA 400 F.O. Otoscope with 1 set (4 pcs.) of reusable tips and 5 of both AllSpec disposable tips 2.5 and 4 mm Ø, spare bulb, **BETA battery handle 2.5V**, hard case [02]

**B-143.10.118**

**Set:** BETA 400 F.O. Otoscope with 1 set (4 pcs.) of reusable tips and 5 of both AllSpec disposable tips 2.5 and 4 mm Ø, **BETA 200 Ophthalmoscope head**, 1 spare bulb each, 1x BETA L rechargeable handle 3.5V, NT300 charger, hard case

**A-153.29.420**

**Set:** BETA 400 F.O. Otoscope with 5 of both AllSpec disposable tips 2.5 and 4 mm Ø, **DELTA 20 Plus Dermatoscope head**, spare bulb for BETA 400, 1x BETA L rechargeable handle 3.5V, NT300 charger, zipper case

**A-160.29.420**



[01]

[02]



**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900  
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077  
Email medical@medsurg.com.au

www.medsurg.com.au  
Fax 07 3859 2995

