

**Mobility and quality have now  
become *one*. The new  
HEINE DELTA*one* Dermatoscope.**



**30**  
YEARS  
1st DERMATOSCOPE



# Small on the outside. Big on the inside.

## The new HEINE DELTAone Dermatoscope.

This mobile masterpiece is taking medical professionals by storm. The crystal-clear image was accomplished using the new achromatic HEINE optical system. What's more, it produces extraordinary colour rendering, due to the quality of our LED<sup>HQ</sup>. Examinations are now particularly comfortable

due to the polarisation allowing doctors to work glare-free and without disrupting reflections. These three features – clarity, LED<sup>HQ</sup> and non-reflective working – are what turn the DELTAone into a vital working tool wherever you go.

Crystal-clear image with the new achromatic HEINE optical system

Fits into any pocket



Glare-free and non-reflective working conditions thanks to polarisation

**LED<sup>HQ</sup>**  
LED NOW IN HEINE QUALITY.

Excellent colour rendering due to LED<sup>HQ</sup>

Digital documentation with the exclusive HEINE DERM app or the HEINE Cube System\*\*

- ∴ **Field of view of real 21mm** (lens 22mm)
- ∴ **Custom light intensity:** 3 degrees of luminous intensity
- ∴ **10x magnification**
- ∴ **Switch between polarised and non-polarised light at the touch of a button** = toggle function. For more information, please visit [heine.com/en/Toggle](http://heine.com/en/Toggle)
- ∴ **No cross-contamination:** autoclavable contact plate with scale

- ∴ **Maintenance-free** due to our dustproof design
- ∴ **Dioptre compensation**
- ∴ **Rechargeable battery with Li-ion technology:** no memory effect; you can safely recharge the battery regardless of its charge state
- ∴ **Wireless connection** via iPhone\* and WLAN with the HEINE Cube System\*\*

**DELTAone Standard Edition** (black-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case

**K-210.28.305**

**DELTAone White Edition** (white-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case

**K-211.28.305**

∴ **Accessories and further articles under [heine.com/en/DELTAone](http://heine.com/en/DELTAone)**

\*The Apple iPhone is not included in the scope of delivery. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. \*\*Not available in all countries.



**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900  
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077  
Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au)

[www.medsurg.com.au](http://www.medsurg.com.au)  
Fax 07 3859 2995

