



Excellence in Human Service  
And Technology

# MINIDOP ES-100VX

POCKET DOPPLER WITH GREAT SOUNDS

## FEATURES

- Pocketable, Hand-held Size
- Great sensitivity
- Wide probe selection for user flexibility  
Clinics: 2, 4, 5, 8 & 10MHz  
Surgical: 8, 10 & 20MHz
- Convenient probe activation button  
turn the unit ON and OFF

## CLINICAL APPLICATIONS

- Blood flow detection



## PROBE SELECTION

### — Standard type —



2MHz Probe



4, 5, 8 & 10 MHz Probe

### — Sterilizable type \*Amplifier required —



Small pencil type



Long pencil type

..and more

## SPECIFICATIONS

<b>Power Supply:</b>	DC 9V, 9V alkaline square battery
<b>Battery life:</b>	Approx. 3 hours
<b>Max speaker output:</b>	300 mW or more
<b>Multi frequency:</b>	2, 4, 5, 8, 10 & 20MHz
<b>Automatic shut-off:</b>	No signal: 3min.
<b>External output:</b>	Headphone Connector: Speaker cut off

<b>Dimensions:</b>	78 (W) x 141 (L) x 27(H) mm (Probe holder not included)
<b>Weight:</b>	240g (Battery and probe are included)
<b>Accessories:</b>	Carrying bag, Gel, Battery
<b>Probe button:</b>	Power ON/OFF
<b>Electrical safety:</b>	Conform to IEC60601-1

*\*Specifications to subject to change.*



Established 1922

## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

