

HEINE® EN 200 Diagnostic Centre

Diagnostic instruments located at your fingertips
in a space-saving design.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



HEINE® EN 200 Diagnostic Centre

Everything in its right place, ready to use: No more tiresome searches for instruments, or manual charging necessary.



Simultaneous use of handles possible due to galvanic separation.

Optimised workflow and energy saving due to Ready to use function.

With the HEINE EN 200 the handle switches on and off automatically when taken off or returned to the unit. Set your brightness level once and it stays there every time you use your instrument.

Exact light conditions and true colour rendering due to continuous brightness control between 100% and 3%.

LED HQ

LED NOW IN HEINE QUALITY.

No unauthorized instrument head removal due to anti-theft mechanism.

Flat, space-saving housing without any openings (no mechanical switch) allows **optimal hygienic reprocessing**.

Connectivity Ready.

The device has an USB interface in order to connect modules such as the HEINE EN 200 BP.



Spiral cord with a length of up to 3.50 m.

Excellent energy efficiency.

With an XHL instrument, typically 58% less power consumption compared to its predecessor EN 100.

Soft Start Technology protects bulbs in the critical switch-on phase (prevents current surge) and prolongs bulb life.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



Whether in a clinic or a practice – the HEINE EN200 can be customised to match your individual needs and can be equipped with a number of HEINE instruments, such as otoscopes, ophthalmoscopes and sphygmomanometers. The wall unit can also be expanded with thermometers of other manufacturers.



BETA F.O. Otoscopes

- ✦ **Solid, all-metal construction.**
Unsurpassed durability.
- ✦ **Exclusive 4,2x or 3x magnification.**
More detail of anatomical structures and easy recognition of foreign bodies.
- ✦ Exclusive continuous brightness control between 100% and 3% (patent pending). **Optimal light intensity is individually adjustable for every examination situation.**
- ✦ **XHL or LED illumination.**



BETA Ophthalmoscopes

- ✦ **Dustproof housing,** maintenance free.
- ✦ Exclusive continuous brightness control between 100% and 3% (patent pending) with one-finger operation. At the same time colour rendering values and homogeneous illumination stay the same. To suit **any examination environment,** and for increased patient comfort **by reducing glare.**
- ✦ **Large, crisp fundus images** with the exclusive Aspherical Optical System from HEINE. Eliminates corneal and iris reflexes.
- ✦ **XHL or LED illumination.**



Accessories

- ✦ **HEINE AllSpec Tip Dispenser.** The right size, visible and instantly available. A clean, new tip for each and every patient.
- ✦ **HEINE EN200 BP Automatic Digital Sphygmomanometer.**
- ✦ **HEINE GAMMA XXL LF.** Analogue Sphygmomanometer, non-mercurial. It is both safe and environmentally-friendly.
- ✦ **HEINE Cuff Basket.**
For storing blood pressure cuffs.

HEINE EN200 BP Automatic Digital Sphygmomanometer with intuitive menu navigation.



HEINE Cuff Basket offers an optimised workflow.



HEINE AllSpec Tip Dispenser. A clean, new tip for each and every patient.



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



EN 200 Wall Transformer

without instrument heads [01]

X-095.12.250**EN 200 Diagnostic Centre**

with EN 200 Wall Transformer, BETA400 LED F.O. Otoscope, BETA200 LED Ophthalmoscope, 1 set of reusable tips (4 pcs.) [02]

A-095.12.208

with EN 200 Wall Transformer, K 180 F.O. Otoscope, K 180 Ophthalmoscope, 1 set of reusable tips (4 pcs.), 50 disposable tips 4mm dia., 50 disposable tips 2.5mm dia. and 1 spare bulb for each instrument

A-095.12.210**AllSpec Tip Dispenser, white**

for AllSpec Disposable Tips [03]

B-000.11.149**EN 200 BP Automatic Digital Sphygmomanometer**

for use with the EN 200 Wall Transformer (with short USB cord), coiled tubing to connect cuffs, without cuffs [04]

M-000.09.201**Cuff Basket**

for cuffs [05]

M-000.09.140**GAMMA XXL LF**

Analogue Sphygmomanometer [06]

M-000.09.323**Wallboard** made of satined acrylic glass for simple wall mounting

for device combination [07]

(50x75x0.8cm, 3.8kg, with pre-drilled holes, mounting material included)

X-095.12.005

for device combination [08]

(50x75x0.8cm, 3.8kg, with pre-drilled holes, mounting material included)

X-095.12.006**EN 200 DIAGNOSTIC CENTRE – CONFIGURATION EXAMPLES:**

[07]



[08]



[01]



[02]



[03]



[04]



[05]



[06]

For additional product configurations and further HEINE Quality Products please visit www.heine.com or give us a call!**HEINE. QUALITY MADE IN GERMANY.**

In 1946 our founder, Helmut A. Heine started HEINE Optotechnik with one goal in mind; to build the highest quality diagnostic instruments available in the world. It was his belief that the highest quality instruments will always best support the health care professional in making the earliest and most accurate diagnosis. For over 70 years, our customers have been able to rely on this promise, and on their HEINE instruments for precise diagnostic performance, unique functionality and long lasting durability.

Today, as a global leader in primary diagnostic instruments with over 500 employees we remain committed to our founding principles. We are still a 100 % family owned and managed company and care deeply about our employees, our customers, and the quality of our products. We continue to develop and manufacture HEINE instruments exclusively at our facilities in Germany, where we combine decades-long experience and craftsmanship with the most modern manufacturing technologies.

HEINE Optotechnik GmbH & Co. KG, Kientalstr. 7, 82211 Herrsching, Germany
Tel. +49 (0) 8152/38-0, Fax +49 (0) 8152/38-202, E-Mail: info@heine.com, www.heine.com



A-000.00.217 e 10/16.
We reserve the right to change specification without notice.
© 2016 HEINE Optotechnik GmbH & Co. KG. All rights reserved.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

